

Reliable versatility to move  
your research forward

**ACT NOW  
FOR  
SAVINGS!**

### Special Package Promotion

Choose the ideal shaker package for your proteomic and genomic application needs and receive a platform and clamps to provide maximized capacity



#### Thermo Scientific™ MaxQ™ Digital Shakers

- Designed to provide reliable results required for bacteria and yeast cultures
- Durable triple eccentric drive backed by 10 year parts warranty
- Maximize capacity of 250 mL and 2 L flasks in a compact footprint
- Digital display with simultaneous viewing of temperature, speed and time with audible and visual alarms








#### Thermo Scientific Max Q refrigerated shakers

- Small scale to large scale units (up to 32kg of loading)
- Cooling from -5°C to -15°C of ambient
- RPM from 15 -500
- Stackable option to increase capacity without the need for additional floor space

**Each package includes shaker, platform and clamps in one part number for easy ordering!**






## Proteomic Packages

- Choose from three shaker configuration options: benchtop, stackable, or large stackable
- Each shaker package is optimized for 250 mL or 2 L flasks to maximize your research
- Incubated units are designed to maintain 37°C making it ideal for bacteria and yeast
- Refrigerated units are designed to hold temperatures at 15°C below ambient to 4°C which is ideal for protein expression

Application package	Shaker description	Accessories description	Pkg Qty	Promo part number	
	MaxQ 4000 with refrigeration	Benchtop Package A	18 x 18" universal platform	1	SHKE40008PA-1CE
			250 mL flask clamps	16	
	MaxQ 6000 with refrigeration	Stackable Package A	18 x 18" universal platform	1	SHKE60007PA-1CE
			2 L flask clamps	6	
	MaxQ 8000 with refrigeration	Stackable Package B	18 x 18" universal platform	1	SHKE60007PB-1CE
			250 mL flask clamps	16	
	MaxQ 8000 with refrigeration	Large Stackable Package A	29.5 x 18" universal platform	1	SHKE80007PA-1CE
			250 mL flask clamps	30	
	MaxQ 8000 with refrigeration	Large Stackable Package B	29.5 x 18" universal platform	1	SHKE80007PB-1CE
			2 L flask clamps	12	

## Genomic Packages

- Each shaker has package configurations optimized for 250 mL or 1 L flasks to maximize your research
- Incubated shaker packages are available in temperatures up to 80°C
- Open air shaker packages are available in a variety of sizes or capacities to fit your lab's needs

Application package	Shaker description	Accessories description	Pkg Qty	Promo part number	
	MaxQ 2000 Digital Shaker	Small Open Air Package A	18 x 18" universal platform	1	SHKE2000GA-1CE
			250 mL flask clamps	16	
	MaxQ 4450 Digital Shaker	Small Open Air Package B	18 x 18" universal platform	1	SHKE2000GB-1CE
			1 L flask clamps	9	
	MaxQ HP420 High Performance Digital Shaker	Small Incubated Package A	13 x 11" universal platform	1	SHKE4450GA-1CE
			1 L flask clamps	4	
	MaxQ HP420 High Performance Digital Shaker	Small Incubated Package B	13 x 11" universal platform	1	SHKE4450GB-1CE
			250 mL flask clamps	9	
	MaxQ HP420 High Performance Digital Shaker	Large Incubated Package A	250/300 mL dedicated platform with 25 pre-installed clamps	1	SHKE420HPGA-1CE
		Large Incubated Package B	1 L dedicated platform with 8 pre-installed clamps	1	

Valid while supplies last. All promotions valid for orders placed between April 15, 2017 and October 15, 2017 by residents of China, India, Singapore, Korea, Thailand, Malaysia, Indonesia, Australia, New Zealand, Taiwan, and Hong Kong. Packages and discounted add-on products must appear on same order to receive discount. "Savings" or "Discount" indicated in any promotion herein represent the maximum percentage difference between the aggregate price for individual purchase of each item contained in the promotional package at the current standard manufacturer's suggested retail price and the manufacturer's suggested retail price for the promotional price for each promotional package. Valid only in connection with purchase of promotional packages identified herein using the part numbers listed. All pricing information and values based on current list prices. All offers subject to Thermo Fisher standard terms and conditions of sale, and Thermo Fisher makes no additional warranties, express or implied, in connection with promotional products. Promotional products will be shipped with purchased products. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. See representative for additional information.

Find out more at [thermofisher.com/shakers](http://thermofisher.com/shakers)

**ThermoFisher**  
SCIENTIFIC