



Autoclave accessories

Odo-Clave® autoclave deodoriser

Neutralise autoclave odors with a choice of four fragrances

Heat-released deodorant pads neutralise autoclave odors with one of four pleasant scents. Simply open foil pouch and place pad in bag prior to autoclaving.

Pads have a three year life span and will last longer if stored in a sealed container.

Item Code	Description	Units	Price
BATH13198.0000	Odo-clave deodorant 32 mm square pad clove fragrance	100pcs	\$197.50
BATH13198.0001	Odo-clave deodorant 32 mm square pad cherry fragrance	100pcs	\$197.50
BATH13198.0002	Odo-clave deodorant 32 mm square pad pine fragrance	100pcs	\$197.50
BATH13198.0003	Odo-clave deodorant 32 mm square pad lemon fragrance	100pcs	\$197.50



Clavies® heat resistant biohazard autoclave/oven gloves

Soft, pliant, all-cotton terry cloth gloves are superb for handling hot objects out of the autoclave or oven; or for other heat protection in the lab, school or industry.

Note: Gloves are not waterproof and should not be used to handle dry ice. CE marked

- Heat-resistant up to 232°C
- Gauntlets in two lengths provide wrist and forearm protection
- Bright orange colour offers high visibility

Item Code	Description	Units	Price
BATH13201.0000	Autoclave glove 33 cm overall length	1 Pair	\$68.00
BATH13201.0001	Autoclave glove 47 cm overall length	1 Pair	\$92.00




Clavies® autoclave safe bag clips

These clips can be opened and closed with one hand to easily and securely close the gathered top of an autoclave bag after it has been filled with wastes.

Stiff plastic design holds the bag neck closed while allowing release of internal pressure that can develop during the autoclave cycle

After autoclaving, the clip is left on the bag for final disposal

Item Code	Description	Units	Price
BATF13190.0100	Autoclave safe bag clips	100pcs	\$35.00



LabServ® autoclave bags

- Designed for maximum protection and convenience
- Steam autoclavable at 121°C
- Made from strong high impact HDPE material

Item Code	Description	Units	Price
LBS3500/502	Biohazard autoclave bags 300 x 560 mm	50pcs	\$25.00
LBS3500/503	Biohazard autoclave bags 410 x 610 mm	50pcs	\$35.00
LBS3500/504	Biohazard autoclave bags 600 x 900 mm	50pcs	\$55.00



Fisherbrand™ white autoclave tape

- Excellent adhesive backing sticks securely
- Removes easily without leaving residue
- White steam-sensitive stripes turn black when passed through a steam sterilisation cycle

Item Code	Description	Units	Price
FSH15-903	Autoclave tape white 55m x 12mm wide	EA	\$15.00
FSH15-904	Autoclave tape white 55m x 19mm wide	EA	\$16.00
FSH15-905	Autoclave tape white 55m x 25.4mm wide	EA	\$18.00



Order Placement:

For customer service, call 1300-735-292
 To fax an order, use 1800-067-639
 To email an order, ordersau@thermofisher.com

Contact Us:

For customer service, call 1300-735-292
 For service and calibration, call 1300-736-767
 To order online: thermofisher.com.au

Find out more at thermofisher.com.au

ThermoFisher
 SCIENTIFIC