



UP TO
**35%
OFF**

Thermo Scientific Reagents Spring Promotion

Parameter	Range	Method	Type	Thermo Orion Catalogue Number	Promotional Price
Aluminium	0.01 - 0.25 mg/L	Erichrome Cyanine R	Powder Pack	AC4P27	\$220.00
Ammonia	1 - 50 mg/L	Salicylate	Tube Test	ACR011	\$147.00
Ammonia	0.02 - 2.50 mg/L	Salicylate	Tube Test	ACR012	\$160.00
Free Chlorine	0.02 - 2.00 mg/L	DPD	Powder Pack	AC4P71	\$32.00
				AC4P71B	\$270.00
Total Chlorine	0.02 - 2.00 mg/L	DPD	Powder Pack	AC4P72	\$32.00
				AC4P72B	\$270.00
COD LR	3 - 150 mg/L	Dichromate	Tube Test	CODL00	\$315.00
COD MR	20 - 1500 mg/L	Dichromate	Tube Test	CODH00	\$315.00
COD HR	200 - 15000 mg/L	Dichromate	Tube Test	CODHP0	\$315.00
Copper	0.05 - 5.00 mg/L	Bicinchoninate	Powder Pack	AC4P29	\$70.00
Iron	0.02 - 3.00 mg/L	Ferro	Powder Pack	AC4P78	\$34.00
Iron	0.02 - 1.80 mg/L	TPTZ	Powder Pack	AC4P79	\$45.00
Manganese	0.01 - 0.70 mg/L	PAN	Powder Pack	AC4P54	\$85.00
Manganese	0.1 - 18.0 mg/L	Periodate oxidation	Powder Pack	AC4P55	\$100.00
Molybdate	0.3 - 40.0 mg/L	Mercaptoacetic Acid	Powder Pack	AC4P42	\$140.00
Nitrite	0.5 - 30.0 mg/L	Chromotropic Acid	Tube Test	ACR007	\$120.00
Nitrate	0.01 - 0.30 mg/L	Diazotization	Powder Pack	AC4P46	\$55.00
Total Nitrogen	0.5 - 25.0 mg/L	Persulphate Digestion	Tube Test	ACD004	\$200.00
Total Nitrogen	5 - 150 mg/L	Persulphate Digestion	Tube Test	ACD007	\$215.00
Oxygen Scavengers (DEHA)	3 - 450 µg/L	Iron Reduction	Liquid	AC4P88	\$165.00
Phosphate	0.06 - 2.50 mg/L	Ascorbic Acid	Powder Pack	AC4P95	\$48.00
Phosphate	0.06 - 5.00 mg/L	Ascorbic Acid	Tube Test	ACR095	\$90.00
Total Phosphorous	0.06 - 3.50 mg/L	Acid Persulphate Digestion	Tube Test	ACD095	\$120.00
Silica	0.01 - 1.60 mg/L	Heteropoly Blue	Powder Pack	AC4P60LR	\$99.00
Silica	1 - 90 mg/L	Silicomolybdate	Powder Pack	AC4P60	\$125.00
Sulphate	5 - 100 mg/L	Turbidimetric	Powder Pack	AC4P82	\$48.00

Terms and Conditions

Prices exclude GST. Offer available until 31st December 2017. Offer not valid with any other discount or offer.

ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-295
Email: InfoWaterAU@thermofisher.com
Visit us online: www.thermofisher.com.au

©2015 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917