



End of Year Sale

Consumables & Equipment Offers

November - December 2017

ThermoFisher
SCIENTIFIC

Contents

ClipTip pipette kits	3	Speedvac concentrator	12
F1 pipette kits	4	Incubators	12
Free serological pipette offer	5	Refrigerated incubators	13
Pipette filler	6	Ovens	13
Disposable filter units	6	Balances	14
Centrifuge tubes	7	Freeze dry systems	14
Desiccator	7	Water activity meters	15
Filter paper	7	Endecotts sieves	15
Gloves	8	Colour meters	16
Syringes	8	Refractometer	17
Digital thermometers	9	Viscometer	17
Dairy testing essentials	9	CO ₂ incubator	18
Centrifuges	10	Refrigerators and freezers	19
Spectrophotometers	10		
Muffle furnaces	11		
Waterbaths	11		

In Australia:

For customer service, call 1300-735-292

To fax an order, use 1800-067-639

To email an order, ordersau@thermofisher.com

Find out more at:

thermofisher.com.au/promotions

ThermoFisher
SCIENTIFIC

All offers subject to Thermo Fisher standard terms and conditions. Errors and omission exempted. Specifications, terms and pricing are subject to change. Please consult your local sales representative for details. © 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Worry-free pipetting is here with the only tip that clips

Thermo Scientific F1-ClipTip GLP kits

Thermo Scientific™ F1-ClipTip™ GLP kits offer a mix of single and multichannel pipettes packaged with compatible tips and a pipette stand. Choose the right kit along with multichannel pipettes for worry-free pipetting and achieve newfound confidence with the Thermo Scientific ClipTip pipetting system.



Key Features & Benefits

- Eliminate the frustration of banging tips on the pipette - tips will lock in place with the lightest touch
- Enable consistent and reproducible pipetting from user to user - higher quality results and more efficient research
- Reduce the risk of loose or leaking tips - once ClipTip pipette tips are on they are sealed and will not loosen until ejected

Item Code	Description	Unit	Price
FNP4701140N	Kit 1: 1x 1-10 µL, 1x 10-100 µL, 1x 100-1000 µL	EA	\$730.00
FNP4701150N	Kit 2: 1x 2-20 µL, 1x 20-200 µL, 1x 100-1000 µL, 1x 30-300 µL 8 channel	EA	\$1,180.00
8 Channel			
FNP4661110N	1-10 µL	EA	\$650.00
FNP4661120N	5-50 µL	EA	\$650.00
FNP4661130N	10-100 µL	EA	\$650.00
FNP4661140N	30-300 µL	EA	\$650.00
12 Channel			
FNP4661150N	1-10 µL	EA	\$720.00
FNP4661160N	5-50 µL	EA	\$720.00
FNP4661170N	10-100 µL	EA	\$720.00
FNP4661180N	30-300 µL	EA	\$720.00

For performance and ergonomic comfort make Thermo Scientific Finnpiquette F1 GLP kits the choice for your lab

Thermo Scientific™ Finnpiquette™ F1 GLP kits offer a variety of single channels in a range of volumes that are packaged with compatible Finntip™ tips and a convenient pipette stand. Choose the right kit required for pipetting in these all-inclusive kits and take advantage of the Finnpiquette F1 performance and lightweight design. Add on a F1 multichannel pipette to your selection to maximise pipetting reliability and comfort.



Key Features & Benefits

- Ergonomically in-touch - light weight handle design combines with low plunger, tip attachment and ejection forces
- Silent effortless volume adjustment and secure volume lock
- Adjustable comfort - 120° adjustable finger rest for both right and left handed operator comfort, allows the hand to relax during pipetting cycles
- Super blow out - 150% increase in air boost to ensure efficient delivery of micro-volumes
- CE marked and 5 year warranty with online registration

Item Code	Description	Unit	Price
FNP4700850N	F1 GLP Kit 1: 1-10 µL micro, 10-100 µL and 100-1000 µL	EA	\$750.00
FNP4700860N	F1 GLP Kit 2: 0.2-2 µL micro, 2-20 µL and 100-1000 µL	EA	\$860.00
FNP4700865N	F1 GLP Kit 3: 10-100 µL, 100-1000 µL and 1-10 mL	EA	\$750.00
FNP4701060N	F1 GLP Kit 4: 2-20 µL, 20-200 µL and 100-1000 µL	EA	\$750.00
8 Channel			
FNP4661000N	1-10 µL micro	EA	\$720.00
FNP4661020N	10-100 µL	EA	\$720.00
FNP4661030N	30-300 µL	EA	\$720.00
12 Channel			
FNP4661040N	1-10 µL micro	EA	\$850.00
FNP4661060N	10-100 µL	EA	\$850.00
FNP4661070N	30-300 µL	EA	\$850.00

Special price and free serological pipettes

Your cells deserve the best

Enhance cell attachment, growth and differentiation with Thermo Scientific™ Nunc™ cell culture treated EasYFlasks™. These flasks receive the Nunclon™ Delta surface treatment that offers maximum adhesion for a broad range of cell types.

- Nunclon Delta certification, including testing on four different cell types, provided in each case
- Two cell lines and one primary cell culture for monolayer formation
- One cell line for cloning efficiency
- Sterilised by gamma irradiation and certified non-pyrogenic
- Ergonomic angled-neck design with a 1/3 turn to open, available in filter or vent/close cap
- Made in Denmark (for the past 60 years!)



Special price and purchase three cases and receive a case of 5 mL or 10 mL new Nunc serological pipettes of your choice free

Quality from start to finish

Reliability • Reproducibility • Scalability • Traceability • Security • Simplicity

Item Code	Description	Unit	Price
Nunclon Delta EasYFlask with filter cap			
NUN156367	25 cm ² EasYFlask with filter cap	Case of 200	\$147.00
NUN156499	75 cm ² EasYFlask with filter cap	Case of 100	\$116.50
NUN159910	175 cm ² EasYFlask with filter cap	Case of 30	\$61.50
Nunclon Delta EasYFlask with vent/close cap			
NUN156340	25 cm ² EasYFlask with vent/close cap	Case of 200	\$120.50
NUN156472	75 cm ² EasYFlask with vent/close cap	Case of 100	\$98.50
NUN159920	175 cm ² EasYFlask with vent/close cap	Case of 30	\$63.00



New Nunc serological pipettes with PET filter

Protect against contamination, non-drip design

Certified sterile, DNase/RNase and human DNA free

Free with every 3 cases of EasYFlask purchased

Item Code	Description	Unit
NUN170355N	Pipette serological 5 mL paper pack sterile	Case of 200
NUN170356N	Pipette serological 10 mL paper pack sterile	Case of 200

Free wing stand

Faster performance with maximum pipetting comfort

Thermo Scientific S1 pipette filler

- Back lit LCD display showing battery life
- Adjustable speed settings
- 2 year warranty with web registration
- Available in 5 colours



Free wing stand MAT9069 with each pipette filler*

*Quote EOY to receive free wing stand

Item Code	Description	Unit	Price
MAT9501	S1 pipette filler white	EA	\$360.00
MAT9511	S1 pipette filler clear	EA	\$360.00
MAT9521	S1 pipette filler blue	EA	\$360.00
MAT9531	S1 pipette filler red	EA	\$360.00
MAT9541	S1 pipette filler green	EA	\$360.00
MAT9057	Hydrophobic filter 0.45 µm sterile for S1 pipettor	Case of 25	POA

Buy 3 and save

Thermo Scientific™ Nalgene™ Rapid-Flow™ sterile disposable filter units with PES membrane provide fast flow rate and high throughput. Low protein binding PES membrane with low extractables are the best choice for sterile filtration of cell culture media, serum, additives and buffers.

Increase sample protection and filtration efficiency

- Filter cup unscrews from storage bottle for easy access to filtrate
- Non-pyrogenic, non-cytotoxic and leakproof
- Gamma radiation sterilised and individually bagged with 5-year shelf life
- 0.2 µm units protect against bacterial contamination



*Price when purchasing 3 cases or more

Item Code	Description	Unit	Price*
NAL564-0020	50 mL, 0.2 µm, 50 mm membrane diameter	12pcs	\$116.00
NAL568-0020	250 mL, 0.2 µm, 50 mm membrane diameter	12pcs	\$132.00
NAL566-0020	500 mL, 0.2 µm, 75 mm membrane diameter	12pcs	\$150.00
NAL567-0020	1 L, 0.2 µm, 90 mm membrane diameter	12pcs	\$206.00

Stock up and save, special offer

Nunc 15 mL and 50 mL polypropylene centrifuge tubes

- DNase and RNase free, USP Class VI, non-pyrogenic, non-cytotoxic and sterile
- Leak proof to protect your samples



Item Code	Description	Unit	1 case	2 cases	3+ cases
NUN339650	15 mL PP conical sterile bulk pack	500pcs	\$95.00	\$85.00	\$75.00
NUN339652	50 mL PP conical sterile bulk pack	500pcs	\$145.00	\$125.00	\$95.00

Break the glass habit and switch to plastic

Thermo Scientific Nalgene desiccator with stopcock

- Clearly view contents with transparent PC cover
- Robust polypropylene body makes it lighter in weight than glass
- Room temperature use only and can be used with or without a vacuum

Item Code	Description	Unit	Price
NAL5310-0250	250 mm with PTFE stopcock, without desiccator plate	EA	\$160.00
NAL5312-0230	230 mm desiccator plate only autoclavable	EA	\$100.00



40% off all Fisherbrand filter papers in stock

Highest quality, competitively priced and reliable, Fisherbrand™ filter paper comes in a wide range including qualitative, quantitative and glass fibre.

Contact our customer service team on 1300 735 292 for more information.



The thinnest chemical resistant synthetic composite disposable glove

The triple layered nitrile and neoprene design of the Microflex® 93-260 glove allows for superior protection against harsh chemicals without sacrificing comfort

- Thin construction provides a high level of tactility in comparison to re-usable gloves
- Extra soft material for greater comfort, flexibility and feel for long periods of use
- Extended cuff for added protection to the wrist and forearm (300 mm total length)
- Powder free and textured fingers for increased grip



Item Code	Description	Unit	Price
ANS93-260-S	Gloves nitrile and neoprene microflex small	50pcs	\$21.00
ANS93-260-M	Gloves nitrile and neoprene microflex medium	50pcs	\$21.00
ANS93-260-L	Gloves nitrile and neoprene microflex large	50pcs	\$21.00
ANS93-260-XL	Gloves nitrile and neoprene microflex x-large	50pcs	\$21.00

The quality syringes you trust - repackaged

Terumo® bulk non-sterile syringes now available

- Less packaging than sterile syringes, reducing unpacking time and waste
- Highly transparent barrels for easy recognition of air bubbles
- Smooth plunger action for greater ease in use
- Latex and natural rubber free



Item Code	Description	Unit
TER6SS80B.01S	Syringe hypodermic leur slip non-sterile 1 mL	4000pcs
TER6SS80B.03L	Syringe hypodermic leur lock non-sterile 3 mL	3600pcs
TER6SS80B.03S	Syringe hypodermic leur slip non-sterile 3 mL	3600pcs
TER6SS80B.05L	Syringe hypodermic leur lock non-sterile 5 mL	2300pcs
TER6SS80B.05S	Syringe hypodermic leur slip non-sterile 5 mL	2300pcs
TER6100410	Syringe hypodermic leur lock non-sterile 10 mL	1250pcs
TER6XS.10ES1	Syringe hypodermic eccentric slip non-sterile 10 mL	1250pcs
TER6SS80B.10S	Syringe hypodermic leur slip non-sterile 10 mL	1250pcs
TER6SS80B.20ES	Syringe hypodermic eccentric slip non-sterile 20 mL	700pcs
TER6SS80B.20L	Syringe hypodermic leur lock non-sterile 20 mL	700pcs

Temperature readings in less than one second

Infrared single and dual laser thermometer guns

- Easy to use - Turn on, point at sample and take temperature reading
- Display backlight permits viewing in low light conditions
- Non-invasive measurement and suitable for solids, semisolids and liquids
- Comes with a Traceable™ calibration certificate valid for 2 years after manufacture



Item Code	Description	Unit	Price
FSH06-664-38	Single laser for maximum temperature readings -60 to 500°C range 0.1°C resolution and accuracy of ± 2% of reading or ± 2°C. 11:1 field of view	EA	\$153.00
FSH06-664-254	Dual laser with Type-K probe for reading comparison -60 to 500°C infrared and Type-K -64 to 1400°C range 0.1°C resolution ± 2% of reading or ± 2°C accuracy 12:1 field of view	EA	\$193.00

For thermometers you can rely on, choose Traceable

Calibrated to NIST standards and NATA recognised

- Suitable for temperature range of -50 to 300°C
- Traceable calibration certificates with every thermometer



Item Code	Description	Unit	Price
FSH14-648-44	Shock and water resistant 0.1°C resolution between -20 to 200°C, otherwise 1°C and ± 1°C accuracy	EA	\$36.00
FSH15-078G	Pocket size 0.1°C resolution between -20 to 200°C, otherwise 1°C and ± 1.5°C accuracy	EA	\$43.00
FSH14-649-100	Waterproof food 0.1°C resolution between -20 to 200°C, otherwise 1°C and ± 1.5°C accuracy	EA	\$29.50

Stock up on dairy testing essentials

Consumables for milk sediment testing



Item Code	Description	Unit	Price
LBSS2000	Milk sediment test funnel	EA	\$330.00
CLDSP20	Filter pad with information card	500pcs	\$250.00
WELLINTINE	Lintine discs	1000pcs	\$195.00

Quick, quiet and convenient

Thermo Scientific Heraeus Pico and Fresco microcentrifuges

Thermo Scientific™ Heraeus™ Pico™ and Fresco™ microcentrifuges combine power, versatility and convenience in a safe, compact, easy-to-use lab instrument. Designed to accelerate your routine sample preparation, these units provide high speed and generous capacity along with outstanding safety and convenience.



- Standard 24 x 1.5/2.0 mL rotor also accommodates mini-preps and spin columns
- Unique ClickSeal® biocontainment rotor lid ensures maximum sample protection
- Intuitive controls and easy-to-read display simplifies operation
- Quiet performance (<56 dBA Pico and <50 dBA Fresco)

Item Code	Description	Unit	Price
THR75002410	Pico 17 ventilated microcentrifuge with 24 x 1.5/2.0 mL rotor & ClickSeal biocontainment lid. Max RCF: 17000 x g	EA	\$2,190.00
THR75002420	Fresco 17 refrigerated microcentrifuge with 24 x 1.5/2.0 mL rotor & ClickSeal biocontainment lid. Max RCF: 17000 x g	EA	\$4,990.00

Spectrophotometers - see the light

Laboratories around the world rely on Thermo Scientific spectrophotometers

- Quick, easy measurements
- Rugged construction
- Reliable results
- Recommended for teaching labs and basic research applications
- Micro-volume options



Item Code	Description	Unit	Price
NDTND-NDLITE	Spectrophotometer NanoDrop Lite for 1 µL samples	EA	\$6,390.00
UNI840-277300	Spectrophotometer GENESYS 30 visible single beam scanning capability / 325-1100 nm / 5 nm bandwidth	EA	\$3,590.00
UNI840-281600	Spectrophotometer Spectronic 200 with onboard software, includes cell and test tube holder	EA	\$2,800.00

Thermo Scientific™ Thermolyne™ premium large muffle furnaces

A rugged furnace with multiple safety features, ideal for industrial applications including ashing and gravimetric analysis

- 1200°C max temperature
- 4 elements in top, bottom and sides for enhanced temperature uniformity
- Built in vent port to remove moisture and contaminants
- Adjustable alarm for over temperature protection
- Door safety switch stops power when door is opened



Item Code	Description	Unit	Price
BARF6020C-33	14 L furnace with single set point controller	EA	\$4,775.00
BARF6020C-33-80	14 L furnace with 8 step programmable controller	EA	\$5,370.00
BARF6030CM-33	14 L furnace with 4 x 16 step programmable controller	EA	\$8,430.00

Thermo Scientific™ Precision™ unstirred waterbaths

Precision redefined for your lab - the ultimate in performance and reliability

- Simple to operate and monitor with icon-based graphical display, and four stored pre-sets
- Small footprint to free up valuable bench space
- Baths supplied with clear polycarbonate lid, diffuser tray, rubber duck and drain hose for tanks 10 L and above
- Adjustable alarm for over temperature protection
- Advanced microprocessor controller is designed for extended functionality
- Temperature range: ambient +5°C to 99.5°C



Item Code	Description	Unit	Price
BARTSGP02	Precision unstirred bath 2 L	EA	\$920.00
BARTSGP05	Precision unstirred bath 5 L	EA	\$930.00
BARTSGP10	Precision unstirred bath 10 L	EA	\$1,095.00
BARTSGP20	Precision unstirred bath 20 L	EA	\$1,260.00
BARTSGP28	Precision unstirred bath 28 L	EA	\$1,340.00

Thermo Scientific™ Savant™ DNA120 Speedvac™ concentrator

Our easy-to-use, small capacity DNA120 is the ideal choice for nucleic acid applications

- Corrosion resistant, PTFE coated chamber provides years of lasting service
- Choice of three drying rates
- Integrated and compact design and extremely quiet operation
- Built-in, oil free vacuum pump provides maintenance-free operation



Item Code	Description	Unit	Price
SAVDNA120-230	DNA 120 Speedvac with 36 x 1.5 mL rotor	EA	\$9,465.00

Thermo Scientific™ Heratherm™ advanced protocol incubators

Whenever you require exceptional performance with superior temperature uniformity and a focus on sample safety, look no further.

- Tight market leading spatial uniformity and stability: No cold spots and allows the whole chamber to be utilised
- Dual convection with fully adjustable fan speed 0-100%
- Broad temperature range from ambient +5 to 105°C
- Ease of use: intuitive user interface, timer and alarms. Stainless steel interior with rounded corners for easy cleaning and corrosion resistant



Item Code	Description	Unit	Price
HEA51028133	IMH60 Advanced protocol incubator 66 L	EA	\$2,850.00
HEA51028134	IMH100 Advanced protocol incubator 104 L	EA	\$3,350.00
HEA51028135	IMH180 Advanced protocol incubator 178 L	EA	\$4,300.00

Thermo Scientific Heratherm refrigerated incubators

Energy efficient, precise temperature control, reliable temperature environment

- Temperature range of +5 to +70°C
- Peltier technology to maintain the precise temperatures for reliable results and excellent temperature uniformity
- Sophisticated timer function with weekly, daily or real time mode and 10 step programmability function
- No hazardous refrigerants and low energy consumption



Item Code	Description	Unit	Price
HEA51031562	IMP180 Refrigerated incubator 178 L	EA	\$7,350.00
HEA51031565	IMP400 Refrigerated incubator 381 L	EA	\$9,350.00

Thermo Scientific Heratherm general protocol ovens

Carry out daily laboratory work with our general protocol ovens, designed with sample protection as top priority



- Choice of gravity convection for gentle drying or mechanical for fast drying and superior chamber spatial uniformity
- Large easy to view vacuum florescent display with on/off timer
- Temperature range up to 250°C, optimally suited for everyday heating and drying applications
- Inner chamber made from corrosion resistant stainless steel with rounded corners for easy cleaning

*15 amp plug required

Item Code	Description	Unit	Price
HEA51028139	OGS60 65 L Gravity convection oven	EA	\$1,800.00
HEA51028140	OGS100 105 L Gravity convection oven*	EA	\$2,400.00
HEA51028141	OGS180 176 L Gravity convection oven*	EA	\$2,900.00
HEA51028151	OMS60 62 L Mechanical convection oven	EA	\$2,760.00
HEA51028152	OMS100 97 L Mechanical convection oven*	EA	\$3,600.00
HEA51028153	OMS180 170 L Mechanical convection oven*	EA	\$4,500.00

Sartorius™ Entris® balances

Get the reliability and quality of a Sartorius balance, but only pay for the features you really need

- **Wide range:** Capacities range from 120 g to over 8 kg with readabilities from 0.1 mg to 1 g
- **Integrated applications:** Including underpan weighing, animal weighing and counting
- **User friendly operation:** Easy navigation with function keys for ease of use, high contrast and backlit display offer excellent visibility
- Optional internal adjustment for optimal weighing accuracy and operating convenience



Item Code	Description	Unit	Price
SARENTRIS124I-1S	Entris balance 120 g x 0.1 mg internal calibration	EA	\$2,400.00
SARENTRIS224I-1S	Entris balance 220 g x 0.1 mg internal calibration	EA	\$2,560.00
SARENTRIS323I-1S	Entris balance 320 g x 1 mg internal calibration	EA	\$1,685.00
SARENTRIS2202I-1S	Entris balance 2200 g x 0.01 g internal calibration	EA	\$1,685.00
SARENTRIS6202I-1S	Entris balance 6200 g x 0.01 g external calibration	EA	\$2,240.00
SARENTRIS5201-1S	Entris balance 5200 g x 0.1 g external calibration	EA	\$1,285.00
SARENTRIS8201-1S	Entris balance 8200 g x 0.1 g external calibration	EA	\$1,375.00

Labconco™ FreeZone™ freeze dry systems

Make anything possible with our broad line of freeze dryers

FreeZone freeze dryers are designed to handle the lyophilisation needs of research and pilot plant laboratories. Collectors available in -50°C and -84°C with capacities ranging from 2.5 to 18 litres. Special 4.5 litre -105°C for samples containing dilute methanol and ethanol.

- Get consistent, high quality results
- Make challenging samples easy
- Plug and play control of drying accessories with many new accessories that are specific for sample types
- All sizes are available with PTFE-coated collectors that can withstand corrosive compounds
- Choice of benchtop or mobile console models

Freeze dryers for today

Built to empower scientists and researchers; your lab will work smarter with the advanced Lyo-Works™ operating system in our latest FreeZone freeze dryers

- On-board collection of sample data with real-time graphing, including collector temperature, vacuum level, and drying accessory parameters
- Alerts and status updates come to you via email
- Exclusive End-Zone™ end point detection accessory notifies you when your samples are dry (no more uncertainty and over drying)
- Keep data at your fingertips. Program parameters and monitor status in real time, and even store your data for repeat protocols and analysis
- Security settings that protect parameters from being changed



Novasina water activity meters

Take the risk out of food manufacturing and reduce the risk of product recalls

Unwanted bacterial or fungal growth poses a health hazard and risk to profit and reputation. Novasina water activity meters can minimise this risk by ensuring your food samples are at an acceptable water activity level. In addition Novasina water activity meters may help compliance with overseas export regulations.



Simple to use

- Load raw sample, close lid, obtain value
- Set “pass or fail” parameters for samples with LabTouch and LabMaster-Neo models
- Fast user calibration procedure with reusable standards
- Exportable results for QC Managers

Wide application capability

- From low to high water content foods
- LabTouch and LabMaster-Neo manages temperature to ensure accurate reporting of sensitive samples

Please contact our customer service team on 1300 735 292 to find out more.

Endecotts sieves

For accurate, dependable results you can rely on Endecotts sieves thanks to their quality of manufacture.

Endecotts sieves are designed and manufactured to offer qualities that make them extremely precise and accurate whilst offering good handling, nesting and strength.

Whether it's a standard test sieve, or something special for such things as coffee or agriculture you'll find all feature meticulous quality in design and manufacture.

Endecotts sieves available from Thermo Fisher Scientific are supplied in a complete range of sizes, depths, materials and certification.

Please contact our customer service team on 1300 735 292 to find out more.



Konica Minolta™ CM-5 spectrophotometer

The CM-5 brings advanced spectrophotometer capabilities in a compact and easy to use format. Ideal for food laboratories, the instrument features a fully guided operation enabling rapid analysis of samples for colour, reflectance, transmittance and whiteness.

Utilise the direct mode for quick and simple colour analysis of solids or the petri dish system for granular or semi solid samples. Transmittance cells can be used to measure liquids in seconds.

Featuring full data analysis and interpretation, on-board method and results memory, fully loaded standards and user management, the CM-5 is the perfect partner for quality control or product development.



CR400/CR410 colour meters

The benchmark portable colour meter where customers require the convenience of a pure L*A*B* or similar colour instrument. Many other colour scales are also available.

Able to do pass/fail, maintain required colour tolerances, store results, use in-situ with various adaptors to suit solids, liquids, pastes and granular material. These are ideal for reliable and fast colour measurements.



CR-20 colour meter

This instrument can be utilised to control colour from raw materials through the process of manufacturing.

It can also quantify the effect that storage or transportation has inflicted on the whiteness of a product.

The CR-20 fits comfortably in the palm of your hand. It is battery-operated and light-weight, so it is extremely portable.



Bellingham + Stanley Refractometer

The use of a benchtop refractometer facilitates convenient and rapid measurement of concentration in a number of liquid and semi-solid samples.

Food applications include the measurement of sugar in grape must and other fresh fruit, processed food, sugar solutions and beverages. Starch, honey and salt solutions may also be measured.

The Bellingham + Stanley RFM340-T Series of high precision Peltier temperature controlled benchtop refractometers, featuring a wide measuring range up to 1.70 RI and new 7" touch-screen user interface, are ideally suited for use in applications within the food and beverage, pharmaceutical, flavours and fragrance industries.

The instrument's new case is ergonomically designed so that any spills are dealt with by the sloping shape of the case or by the PEEK spill barrier that surrounds the shallow, easy-clean stainless steel prism dish. The prism dish is shaped to hold sufficient sample for measurement whilst regaining a large enough surface area for rapid temperature stabilisation to occur.



Please contact our customer service team on 1300 735 292 to find out more.

Fungilab viscometers for quality control

Directly compatible with Brookfield methods

Easy to use and supplied complete with

- Temperature probe
- Quick PPR spindle connect-disconnect system
- Audible warnings if unit is being operated outside acceptable parameters
- USB connectivity and data export software
- 2 year warranty
- Automatic instrument self diagnosis upon start up

Please contact our customer service team on 1300 735 292 to find out more.



Double your capacity and save

The Thermo Scientific™ HERAcell™ VIOS CO₂ incubator is specifically designed to deliver the performance, reliability, ease of operation, and value required to support a range of culturing needs, from basic research to demanding, leading-edge applications so you're ready for whatever comes next.

THRIVE active airflow

- In-chamber fan gently distributes clean, humidified air throughout the chamber ensuring homogeneous conditions
- Fast recovery of all parameters - ten minutes or less for stable culturing conditions
- Auto-stop function disables fan operation during door openings to minimise air exchange
- In-chamber HEPA continuously filters the entire chamber air volume every 60 seconds

Steri-Run, high temperature disinfection cycle

- Overnight cycle simplifies cleaning and eliminates the need for separate autoclaving of parts
- 180°C cycle assures uniform, 6-log sterilisation on all chamber surfaces eliminating biological contaminants

In situ sensors and probes

- Dual temperature probes for over temperature protection with operational back up
- Choice of thermal conductivity (TC) with rH compensation or single beam, dual wavelength infrared (IR) CO₂ sensor



Save now

Our double capacity incubator package including 2 x HERAcell VIOS 160i CO₂ incubators with stainless steel interior and TC sensor and free Stacking adaptor with roller dolly

* An additional charge of \$500 (excl GST) applies per purchase order for sensitive freight delivery

Item Code	Description	Unit	Price*
HEA50145502	HERAcell VIOS 160i CO ₂ Incubators with stainless steel interior and TC180 CO ₂ Sensor Package includes: 2x CO ₂ incubators, 1x stacking adaptor and 1x roller dolly	EA	\$15,990.00

Protect routine samples:

Selecting the right cold storage equipment has become a critical decision for today's laboratory

Laboratories today are recognising the key differences between standard refrigerators and freezers and those designed specifically for protecting precious samples.

There are three main types of refrigerators and freezers for three specific needs:

Those built for laboratory applications

These refrigerators and freezers are designed to store high value, critical samples, reagents, enzymes and cells. They're specifically engineered to meet the demanding temperature requirements of these valuable items. Sample integrity and protection is paramount.

Those built for domestic use and those built for retail and catering sector

These products are designed for household (such as eggs, meat and cheese) and commercial groceries (such as bottled water, and bulk storage food). If the units' temperature fluctuates outside the required limits, these items will be inexpensive to replace.

The right choice for sample storage

The Thermo Scientific GP Series includes medium and large capacity refrigerators and freezers designed for routine sample protection. With capacity ranging from space-saving to spacious, our new GP Series lab refrigerators and freezers meet the demands of today's laboratory requirements. All models contain sample protection features designed for the rigorous demands of the laboratory environment.

All GP Series models deliver the sample protection, performance, security and quality you have come to depend on from us:

- Integrated controller
- Digital display
- Low energy consumption
- Choice of glass or solid door (refrigerators only)
- Forced air circulation
- Automatic defrost
- HCFC- and CFC-free insulation
- Reversible doors (R/F 400 and R/F 700 models only)
- High and low temperature alarms
- Door ajar alarm
- Remote alarm contacts
- Standard door locks
- Standard casters
- Access ports
- Interior fluorescent lighting



 Request a quotation at:
thermofisher.com.au/gpscoldstorage



End of Year Clearance

Consumables & Equipment Offers

November - December 2017

ThermoFisher
SCIENTIFIC

Pipettes

Item Code	Description	Unit	Stock on Hand	Price
EPP3120000.020T/I	Research Plus Pipette 0.5-10 µL	EA	19	\$270.00
EPP3120000.038T/I	Research Plus Pipette 2-20 µL	EA	3	\$270.00
EPP4961000.012T/I	Varispenser plus 0.5-2.5 mL	EA	1	\$270.00
EPP4910000.077	Reference pipettor variable 500-2500 µL	EA	2	\$400.00
FNP4641070	Finnpipette F1 10-100 µL	EA	1	\$280.00
FNP4641100	Finnpipette F1 100-1000 µL	EA	3	\$280.00
FNP4641130	Finnpipette F1 5-50 µL micro	EA	1	\$280.00
FNP4641170	Finnpipette F1 ClipTip 1-10 µL	EA	8	\$280.00
FNP4641190	Finnpipette F1 ClipTip 5-50 µL	EA	1	\$280.00
FNP4641200	Finnpipette F1 ClipTip 10-100 µL	EA	3	\$280.00
FNP4641210	Finnpipette F1 ClipTip 20-200 µL	EA	3	\$280.00
FNP4641220	Finnpipette F1 ClipTip 30-300 µL	EA	1	\$280.00
FNP4651080	Finnpipette F1 fixed 1000 µL	EA	2	\$170.00
FNP4651110	Finnpipette F1 fixed 5000 µL	EA	1	\$170.00
FNP4661020	Finnpipette F1 8-channel 10-100 µL	EA	2	\$560.00
FNP4661040	Finnpipette F1 12-channel 1-10 µL micro	EA	1	\$840.00
FNP4661060	Finnpipette F1 12-channel 10-100 µL	EA	1	\$820.00
FNP4661070	Finnpipette F1 12-channel 30-300 µL	EA	3	\$670.00
FNP4661120	Finnpipette F1 ClipTip 8-channel 5-50 µL	EA	8	\$520.00
FNP4661130	Finnpipette F1 ClipTip 8-channel 10-100 µL	EA	2	\$520.00
FNP4661140	Finnpipette F1 ClipTip 8-channel 30-300 µL	EA	11	\$520.00
FNP4661150	Finnpipette F1 ClipTip 12-channel 1-10 µL	EA	1	\$625.00
FNP4661160	Finnpipette F1 ClipTip 12-channel 5-50 µL	EA	4	\$625.00
FNP4661170	Finnpipette F1 ClipTip 12-channel 10-100 µL	EA	2	\$625.00
FNP4661180	Finnpipette F1 ClipTip 12-channel 30-300 µL	EA	2	\$625.00
FNP4670020	E1 ClipTip electronic pipettor 2-125 µL	EA	3	\$895.00
FNP4670040	E1 ClipTip electronic pipettor 15-1250 µL	EA	1	\$895.00
FNP4671000	E1 ClipTip electronic 8-channel 0.5-12.5 µL	EA	2	\$895.00
FNP4671040	E1 ClipTip electronic 8-channel 2-125 µL	EA	1	\$1,420.00
FNP4671070	E1 ClipTip electronic 8-channel 10-300 µL	EA	4	\$1,420.00
FNP4671080	E1 ClipTip electronic 12-channel 10-300 µL	EA	1	\$2,040.00
FNP4671100	E1 ClipTip electronic 8-channel 15-1250 µL	EA	1	\$895.00
FNP4672010	E1 ClipTip electronic 8-channel 0.5-12.5 µL with adjustable tip spacing 4.5 to 14.2 mm	EA	1	\$1,260.00
FNP4672050	E1 ClipTip electronic 8-channel 2-125 µL with adjustable tip spacing 9 to 14.2 mm	EA	1	\$1,570.00

FN4672080	E1 ClipTip electronic 8-channel 10-300 µL with adjustable tip spacing 9 to 14.2 mm	EA	1	\$2,200.00
FN4672100	E1 ClipTip electronic 8-channel 50-1250 µL with adjustable tip spacing 9 to 14.2 mm	EA	3	\$2,000.00

Pipette tips

Item Code	Description	Unit	Stock on Hand	Price
FSB02-707-406	Pipette tips 100-1250 µL blue racked	960pcs	19	\$29.00
FSB02-707-410	Pipette tips 5-300 µL sterile racked	960pcs	15	\$36.00
FSB02-707-411	Reload tips 5-300 µL non sterile	960pcs	99	\$29.00
FSB02-707-431	Filter tips 10-100 µL sterile bevelled racked	960pcs	43	\$45.00
FSB02-707-432	Filter tips 2-20 µL sterile bevelled racked	960pcs	61	\$45.00
FSB02-707-442	Filter tips 0.1-10 µL sterile extend length racked	960pcs	24	\$50.00
EPP0030073.282	Pipette tips epTIPS Set 20-300 µL rack and refill	480pcs	6	\$56.00
EPP0030073.347	Pipette tips epTIPS Set 2.5 mL reload	240pcs	7	\$56.00
EPP0030073.827	Pipette tips epTIPS Reloads 20-300 µL	960pcs	12	\$80.00
EPP0030075.048	Tips pipette epTIPS Racks 20-300 µL	480pcs	59	\$65.00
QSP117	Pipette tips 5-250 µL yellow bulk	1000pcs	29	\$30.00

Plastic tubes and containers

Item Code	Description	Unit	Stock on Hand	Price
LBSUC2005X	Container 30 mL PS flat base label yellow	500pcs	34	\$165.00
LBSUC2006N	Container 30 mL PS natural PE cap separate	800pcs	16	\$85.00
LBS32250R	Container 250 mL PS round red PE cap	140pcs	10	\$55.00
LBSHEC110AQ	Cassette embedding aqua	4000pcs	4	\$115.00
LBSBC16G	Cap overcap green for 16 mm blood tubes	1000pcs	102	\$18.00
LBS9615CAPN	Cap PE natural for 9615 tube	2000pcs	1	\$140.00
LBSP20001X	Tube only PS round base 16 x 150 mm	800pcs	2	\$70.00
LBSSP5201	Tube 2 mL PP	2000pcs	8	\$65.00
LBSSP5204	Tube 2 mL PP label white PE screw cap	2000pcs	3	\$255.00
LBSSP2001	Tube 5 mL PS flat base	2000pcs	2	\$66.15
LBS514N	Tube 5 mL PP 11 x 75 mm screw top	4000pcs	7	\$120.00
LBSCF0001	Tube 10 mL PS V base plug tube	1000pcs	18	\$70.00
LBSSCF001	Tube only 10 mL PS V base	1000pcs	50	\$65.25
LBSP00006	Tube 15 mL PS natural PE plugcap separate	2000pcs	7	\$325.00
COS4CC	Container 4 L urine HDPE	36pcs	21	\$46.00
COSA20W	Jar white HDPE 700 mL 83 mm wide neck	140pcs	1260	\$133.00
COSBPS20	Pail 20 L white tamper evident super baroda	EA	12	\$6.30

Beakers and bottles

Item Code	Description	Unit	Stock on Hand	Price
BOKBOE5011624	Beaker 100 mL glass tall form graduated	10pcs	1673	\$9.50
BOKBOE5010668	Beaker 3 L glass low form graduated	EA	14	\$15.45
LBS20202-1L	Beaker glass 1 L Berzelius tall form double scale	EA	70	\$3.30
KIM14000-4000	Beaker glass 4000 mL low form	EA	27	\$64.40
LBS20202-400	Beaker glass 400 mL Berzelius tall form double	EA	750	\$1.70
PYR1015/22D	Beaker glass 1000 mL tall form graduated with spout	EA	78	\$7.50
KIM14000-10	Beaker glass 10 mL low form	EA	40	\$6.15
KIM14000-20	Beaker glass 20 mL low form	EA	42	\$5.25
PYR1020/06D	Beaker glass 250 mL conical Phillips	EA	20	\$7.70
KIM14030-300	Beaker glass 300 mL tall form spouted	EA	76	\$4.10
AZLBRE006	Bottle density 100 mL adjusted Gay Lussac type	EA	10	\$89.80
SCH21801915	Bottle laboratory 20 L graduated with PP cap	EA	3	\$676.50
BOKBOE5080865	Bottle laboratory 10 L glass blue PP screw cap	EA	2	\$104.55
WHE219371	Bottles diagnostic amber 10 mL	480pcs	1	\$850.65
WHE219362	Bottles diagnostic screw thread 20 mL	480pcs	1	\$846.95
LOM1401/500P	500 mL reagent bottle clear	EA	186	\$2.00

Flasks

Item Code	Description	Unit	Stock on Hand	Price
QFTFE100/1	Flask erlenmeyer glass 100 mL 14/23 socket	EA	50	\$11.50
QFTFE100/2	Flask erlenmeyer glass 100 mL 19/26 socket	EA	48	\$11.40
PYR3151/16	Flask volumetric 200 mL class A 14/23 PE stopper	EA	27	\$12.30
LBS20401-02-1L	Flask erlenmeyer 1000 mL narrow neck graduated	EA	83	\$5.70
KIM26650-2000	Flask erlenmeyer 2000 mL wide neck	EA	46	\$16.30
GLS071.202.18	Flask erlenmeyer 250 mL glass with 19/26 socket	EA	956	\$5.20
KIMKC26720-250	Flask erlenmeyer 250 mL plastic coated wide neck	6pcs	4	\$249.30
LBS20401-02-3L	Flask erlenmeyer 3000 mL narrow neck graduated	EA	42	\$15.45
KIMKC26720-500	Flask erlenmeyer 500 mL plastic coated wide neck	EA	5	\$215.95
QFTFE250/2	Flask erlenmeyer glass 250 mL 19/26 socket	EA	50	\$16.35
QFTFE500/3	Flask erlenmeyer glass 500 mL 24/29 socket	EA	38	\$19.10
BOKBOE5120154	Flask filter 1 L glass with glass hose connection	EA	7	\$26.85
LBS40620-250	Flask volumetric 250 mL PE stopper Class A	EA	609	\$3.30
LBS40620-25	Flask volumetric 25 mL PE stopper Class A	EA	75	\$3.82
GLS129.202.10	Flask volumetric glass Class A 2 L 29/32 stopper	EA	27	\$14.05
LOM1621/25G	Volumetric Flask 25 mL with glass stopper	EA	2	\$2.00

Funnels

Item Code	Description	Unit	Stock on Hand	Price
SCH2135138	Funnel 70 mm short stem	10pcs	130	\$193.00
FSBFB68711	Funnel filter 125 mL porosity 4 frit glass	EA	10	\$62.45
FSBFB68714	Funnel filter 500 mL porosity 3 frit glass	EA	10	\$118.80
GLS169.202.04	Funnel pressure equalising with 24/29 socket 250 mL	EA	6	\$19.50
FSBFB51759	Funnel separating 1000 mL pear PTFE stopcock	EA	8	\$77.70
PYRFST321	Funnel separating 100 mL conical PTFE key 2 mm bore	EA	5	\$54.70
PYRFST374	Funnel separating 2000 mL conical PTFE key 6 mm bore	EA	2	\$134.30
PYRFST331	Funnel separating 250 mL conical PTFE key 3 mm bore	EA	3	\$59.50
PYRFST342	Funnel separating 500 mL conical PTFE key 4 mm bore	EA	5	\$85.10
PYRFST311	Funnel separating 50 mL conical PTFE key 2 mm bore	EA	5	\$51.00
PYRFST221	Funnel separating pear shape 100 mL 19/26	EA	4	\$47.00
PYRFST242	Funnel separating pear shape 500 mL 24/29	EA	4	\$69.80
PYRFST252	Funnel separating pearshape 1000 mL 24/29	EA	4	\$76.85

Pipettes and cylinders

Item Code	Description	Unit	Stock on Hand	Price
BOKLMS63336165	Pipette measuring glass Class AS 2 mL	EA	6	\$1.75
GLS124.202.01	Pipette measuring glass graduated Class AS 1 mL	EA	66	\$1.00
GLS124.202.02	Pipette measuring glass graduated Class AS 2 mL	EA	12	\$1.20
GLS124.202.03	Pipette measuring glass graduated Class AS 5 mL	EA	80	\$1.10
KIM72100-1110	Pipette serological 1 mL x 0.1 glass	1000pcs	48	\$338.80
FSBFB50204	Pipette volumetric glass class A 100 mL	2pcs	98	\$27.00
FSBFB50202	Pipette volumetric glass class A 50 mL	2pcs	132	\$15.60
GLS122.202.08	Pipette volumetric glass Class AS 100 mL	EA	15	\$2.80
GLS122.202.04	Pipette volumetric glass Class AS 10 mL	EA	256	\$1.30
GLS122.202.06	Pipette volumetric glass Class AS 25 mL	EA	34	\$1.40
GRF1.100-37-404	Pipette graduated 10 mL x 0.1 Class B	EA	116	\$2.60
GRF1.100-32-404	Pipette graduated 5 mL x 0.05 Class B	EA	138	\$2.50
VOLJ9000PV/250	9 Pasteur Pipette Cotton	250pcs	3	\$25.20
GRF1.350-37-02	Cylinder measuring 10 mL with PP base Class B	EA	143	\$1.80
PYR3201/12	Cylinder measuring 250 mL glass class A	EA	36	\$17.50
GLS138.202.01A	Cylinder measuring 10 mL glass hexagonal base	EA	136	\$2.40
GRF1.230-51-02	Pipette volumetric 100 mL 1 mark Class AS	EA	1	\$9.60

Glass jars and vials

Item Code	Description	Unit	Stock on Hand	Price
COSH15131/P	Jar 250 mL Clear Glass Round pallet	2400pcs	9	\$977.95
FSB02-911-460	Jar clear graduated 125 mL with PTFE line cap	24pcs	3	\$49.50
FSB02-911-458	Jar clear graduated 30 mL with PTFE lined cap	48pcs	3	\$58.00
FSB02-911-459	Jar clear graduated 60 mL with PTFE lined cap	48pcs	3	\$63.20
WHE986492	Vial sampule 6 mL with urea cap	200pcs	49	\$83.70
WHE986540	Vial scintillation glass 20 mL 22 mm PP cap	500pcs	25	\$224.00

Test tubes and tissue grinders

Item Code	Description	Unit	Stock on Hand	Price
GLS082.202.07	Tube culture 60 mL 25 x 200 mm with screw cap	50pcs	2250	\$55.00
SAMG050/27	Tube specimen soda glass PE stopper 75 x 19 mm	100pcs	102	\$29.20
WHE357574	Tissue grinder dounce 15 mL stainless	EA	1	\$626.20
WHE357546	Tissue grinder dounce 40 mL glass	EA	2	\$341.10
WHE357576	Tissue grinder dounce 40 mL stainless	EA	1	\$786.80
WHE357422	Tissue grinder tenbroeck 2 mL	2pcs	12	\$191.60
SCH60101.01	Centrifuge tube glass 15 mL	100pcs	10	\$183.80

Filters and filter paper

Item Code	Description	Unit	Stock on Hand	Price
WHA2800-266	Extraction thimble ID 26 x L 60 mm	25pcs	17	\$63.20
WHA1882-047	Filter circles EPM 2000 4.7 cm	100pcs	80	\$36.65
WHE635233	Filter flask 2000 mL No. 8 stopper joint	EA	5	\$227.80
FLT1926-240F	Filter paper 240 mm No.1926 qualitative	100pcs	14	\$70.30
WHA1441-025	Filter Paper Grade 41 Circles 25 mm	400pcs	2	\$65.70
WHA1542-185	Filter paper grade 542 18.5 cm	100pcs	13	\$67.00
WHA10347577	Filter paper grade sharkskin 29 cm	100pcs	18	\$43.30
HOLCREPED203MM	Paper filter creped 203 mm	100pcs	15	\$18.00
HOLTHNWHITE203MM	Paper filter thin white 203 mm	100pcs	94	\$21.00
PACAP45641	Syringe filter Acrodisc Psf	50pcs	1	\$688.20
MILUFC9050596	Amicon Ultra 15 50 kDA 15 mL	96pcs	2	\$1,224.70
MILSSWP01300	Membrane filter 3.0 µm 13 mm plain MCE	100pcs	21	\$134.70

Crucibles

Item Code	Description	Unit	Stock on Hand	Price
LOMC1	Crucible porcelain high wall with lid 15 mL	45pcs	180	\$35.00
BHAS1460	Crucible silica 80 mL low form translucent	EA	9	\$33.20
BHAS1470	Crucible silica 150 mL low form translucent	EA	14	\$35.35
EXP79MF/7	Crucible medium form porcelain 15 mL 35 mm diameter	10pcs	7	\$100.60

Gloves

Item Code	Description	Unit	Stock on Hand	Price
LOMGNBPF1000MT	Gloves x-small nitrile powder free	Case of 10	4	\$93.10
SUPGLLPF1003MT	Gloves large latex powder free	100pcs	65	\$7.65
S-M10115123	Gloves large latex powder free	100pcs	52	\$8.50

Storage boxes

Item Code	Description	Unit	Stock on Hand	Price
LBSSLD054	Slide storage box and lid for 450 slides	EA	161	\$18.90

Racks

Item Code	Description	Unit	Stock on Hand	Price
SIE135	Test tube rack PP 8 place for 20 mm diameter tubes	EA	25	\$10.70
BATF18745.0002	Test tube rack PP 40 place for 20 mm diameter tubes	EA	137	\$37.35
SIE566/03	Universal tube rack PP 40 place 20 mm holes	EA	18	\$17.40

Tubing

Item Code	Description	Unit	Stock on Hand	Price
TYG116-0549-070	Tygon pump tubing 0.76 mm ID 455 mm length	10pcs	21	\$35.35
TYG116-0549-080	Tygon pump tubing 0.89 mm ID 455 mm length	10pcs	56	\$34.60
TYG116-0549-100	Tygon pump tubing 1.14 mm ID 455 mm length	10pcs	74	\$32.90
TYGAAC00007	Tygon R-3603 1/8 in ID 1/4 in OD 50 ft/roll	EA	7	\$25.00

Tubes

Item Code	Description	Unit	Stock on Hand	Price
SIE88323	Test tube 10 mL PP 16 x 100 mm with rim	2000pcs	12	\$97.60
SIE1298	Micro tube 0.5 mL PP Eppendorf type tube	1000pcs	72	\$19.75

Measuring jugs and cylinders

Item Code	Description	Unit	Stock on Hand	Price
SIELVL646/95	Cylinder measuring 10 mL PMP transparent tall form	EA	319	\$2.25
SIE1160	Jug measuring 5 L PP graduated short form	EA	12	\$9.15

Thermometers

Item Code	Description	Unit	Stock on Hand	Price
BOKBOE305	Thermometer electronic wide range -50 to 300°C	EA	15	\$34.50
WIK538.12276	Magnetic surface hotplate thermometer 10 to 400°C	EA	5	\$125.65
TNOGE5045-01	Dial thermometer 0 to 250°C S/S 130 mm length	EA	8	\$27.30

Plastic bottles

Item Code	Description	Unit	Stock on Hand	Price
BKH0250SL38BWOODBW200	Bottle 500 mL HDPE natural 38 mm snaplock neck	220pcs	6	\$58.80
BKH1000SL38BW120	Bottle 1 L square HDPE natural 38 mm snaplock neck	120pcs	14	\$39.00

Basins

Item Code	Description	Unit	Stock on Hand	Price
LOMC-50	Basin porcelain 75 mL round bottom 82 mm diameter	EA	99	\$0.70
LOMC56	Basin porcelain 180 mL round bottom with spout	EA	61	\$1.15
LOMC57	Basin porcelain 200 mL round bottom with spout	EA	44	\$1.15

Stoppers

Item Code	Description	Unit	Stock on Hand	Price
LBSM171C	Rubber stopper solid 10 mm bottom 14 mm top diameter	EA	4636	\$0.40
LBSM180	Rubber stopper solid 19 mm bottom 22 mm top diameter	EA	2440	\$0.55
LBSM182	Rubber stopper solid 22 mm bottom 25 mm top diameter	EA	696	\$0.70
LBSM798A	Rubber stopper solid 6 mm bottom 10 mm top diameter	10pcs	200	\$3.40

Pipette tips

Item Code	Description	Unit	Stock on Hand	Price
SIE956	Tip micro pipette 0.5-10 µL gilson type	100pcs	5	\$13.70
SIE960	Tips 2-200 µL PP universal yellow	1000pcs	62	\$7.75
SIE964	Tips 1000-5000 µL PP neutral	100pcs	9	\$9.45

Miscellaneous laboratory products

Item Code	Description	Unit	Stock on Hand	Price
SIE418	Valve non return PP 8-10 mm nominal OD	10pcs	22	\$13.95
SIE455	Watch glass PP 99 mm diameter	EA	26	\$3.05
SIE529	Y connector 120 degree PP 8-10 mm OD tapered bore	EA	26	\$2.10
BATF37886.0000	Pipette aid vikem vinyl blue bulb	6pcs	93	\$91.70
DRA8101991	Tube hydrogen sulphide 0.2/b 0.2 - 6 ppm	10pcs	5	\$156.00
LOMSYR005	Syringe 5 mL and plunger without rubber gasket	1000pcs	24	\$38.30
SIFC154	Mortar and pestle porcelain cc154 75 mL 60 mm	EA	52	\$0.90

Serological pipettes

Item Code	Description	Unit	Stock on Hand	Price
LBSCN857110	Pipette serological 10 mL PS individually wrapped	400pcs	234	\$105.00
SLN40101K	Pipette serological 1 mL individually wrapped plugged	Case of 1000	14	\$276.00
NAL6803-0001	Pipette rack with stand acrylic	EA	6	\$132.00

Cell cultureware

Item Code	Description	Unit	Stock on Hand	Price
NUN153066-PROMO	Culture dish with cell culture surface 35 x 10 mm	Case of 500	13	\$51.50
NUN267313	96 well Edge plate nontreated sterile	Case of 50	3	\$292.00
NUN174904	UpCell culture dish 3.5 cm with lid & membrane	Case of 30	2	\$705.00
NUN177380	Lab-Tek chamber slide 2 well	Pack of 16	11	\$115.00
NUN177399	Lab-Tek chamber slide 4 well	Pack of 16	19	\$125.00
NUN155411	Lab-Tek chambered coverglass 8 well	Pack of 16	54	\$122.00
NUN170920	SlideFlask 9 cm ² cell culture surface sterile	Case of 50	5	\$518.00
NUN164588	96 well plate PS CVG cell culture black	Case of 30	1	\$861.50
NUN164327	Cell Factory 10 tray PS nunclon sterile	Case of 2	4	\$252.50
NUN173239	Cell factory active gassing 4 chamber	Case of 10	2	\$956.50
NUN173238	Cell factory active gassing 10 chamber	Case of 6	2	\$1,055.50

Cryotubes and cryoboxes

Item Code	Description	Unit	Stock on Hand	Price
NUN374220	Cryobank tubes 5.0 mL 2D coded racked sterile	Case of 480	3	\$487.50
MAT4247	1.4 mL PP RB tube Latch Rack of 96 tubes/case	Case of 10	7	\$185.00
NAL5025-0505	Cryobox PC 5 x 5 array for 25 1.2 or 2 mL vials	Pack of 8	95	\$58.00
NAL5050-0001	Box storage PC with clear cover 133 x 133 x 51 mm	Pack of 4	46	\$42.00

Genomic and immuno plates

Item Code	Description	Unit	Stock on Hand	Price
ABGAB-0661	Storage plate natural 2.2 mL 96-well 50pcs	EA	35	\$245.00
NUN6309	Module plate 2HB F12	Case of 100	97	\$553.75

Bioproduction vials and bottles

Item Code	Description	Unit	Stock on Hand	Price
NAL312158-0021	Closure 20-415 HDPE 20 mm for serum vial white	Case of 2000	27	\$542.50
NAL322032-0010	Serum vials 10 mL continuous thread PETG	Case of 1260	4	\$395.00
NAL322032-0020	Serum vials 20 mL continuous thread PETG	Case of 612	38	\$225.00
NAL382019-0250	Bottle 250 mL certified clean square media PETG	Case of 48	9	\$362.50

NAL382019-1000	Bottle 1000 mL certified clean square media PETG	Case of 24	5	\$361.25
NAL2015-1000	Bottle 1 L PC square narrow mouth with PP screw cap	Case of 24	13	\$350.00

Jerricans

Item Code	Description	Unit	Stock on Hand	Price
NAL2240-0015	Jerrican 6 L HDPE with integral spout and handle	EA	25	\$33.00
NAL2240-0025	Jerrican 10 L HDPE with integral spout and handle	Case of 6	6	\$268.00

Centrifuge tubes and bottles

Item Code	Description	Unit	Stock on Hand	Price
NUN362696	EZ flip PP 50 mL conical centrifuge tube bulk pack	Case of 500	34	\$189.00
NAL3105-0015	Tube centrifuge 15 mL conical PC	Case of 10	72	\$23.50
NAL3114-0050	Centrifuge tube oak ridge 50 mL FEP	Pack of 2	13	\$94.00
NAL3115-0030	Centrifuge tube oak ridge PSF 30 mL	Pack of 10	8	\$72.00
NAL3115-0050	Centrifuge tube oak ridge PSF 50 mL	Pack of 10	9	\$100.00
NAL3117-0120	Centrifuge tube cylindrical 12 mL PC	Pack of 10	17	\$30.00
NAL3117-0160	Centrifuge tube round PC 16 mL	Pack of 10	24	\$21.00
NAL3127-0250	Centrifuge bottle 250 mL Teflon FEP	Case of 6	11	\$825.00
NAL3138-0050	PC tube 50 mL centrifuge oak ridge with seal	Pack of 10	35	\$135.00
NAL3139-0030	Tube centrifuge PPCO 30 mL oak ridge	Pack of 10	8	\$97.00
NAL3140-1006	Centrifuge bottle PC Super speed 1 L SCA	Case of 6	3	\$685.00

Carboys

Item Code	Description	Unit	Stock on Hand	Price
NAL2211-0050	Carboy 20 L HDPE rectangular with handle and PP cap	EA	99	\$99.00
NAL2234-0020	Carboy 10 L LDPE wide mouth with handles	EA	20	\$58.00
NAL2234-0050	Carboy 20 L LDPE wide mouth with handles	Case of 4	1	\$415.00
NAL2235-0020	Carboy 10 L PP wide mouth with handles autoclavable	EA	45	\$75.00
NAL2319-0130	Carboy 50 L PP w PP spigot & screw cap autoclavable	EA	4	\$330.00
NAL2323-0015	Carboy lowboy 15 L HDPE with PP spigot	EA	9	\$92.00
NAL2324-0015	Carboy lowboy 15 L PP with PP spigot autoclavable	Case of 4	5	\$550.00
NALDS2227-0020	Stirrer carboy magnetic for 10 and 20 L	EA	7	\$305.00

General labware

Item Code	Description	Unit	Stock on Hand	Price
NAL6565-0001	Carrier bottle 4 in 1 EZ tote blue versatile	Case of 4	10	\$158.00
NALDS0390-0070	Filling bell PC 70 mm with PP tubing adaptor	EA	47	\$52.00
NAL2752-9125	Bottle dropper 4 mL LDPE with assorted coloured	Case of 25	62	\$32.50
NAL4301-0125	Funnel separating 125 mL Teflon autoclavable	EA	32	\$115.00

General lab equipment

Item Code	Description	Unit	Stock on Hand	Price
BAR4637-1CEQ	Thermal rocker incubator temp ambient +/-5 to 70°C	EA	1	\$2,395.00
BAR50087930	Magnetic stirring bar round 12 mm D x 30 mm L	5pcs	5	\$85.55
BAR88880030	Heating and cooling block heater dry bath	EA	1	\$1,869.00
BAR88880150	Temperature probe PT100 A class 400°C	EA	1	\$256.00
BARHP131720-33Q	Super-Nuova hotplate ceramic top 18.4 x 18.4 cm	EA	1	\$480.00
ELEEM5000/CE	Mantle heating 5L controlled and earthed	EA	1	\$775.00
ELEEMX5000/SCE	Mantle heating 5000 mL spillproof	EA	1	\$956.00
IKA2983100	A 11.5 Spare grinding chamber (80 mL)	EA	1	\$275.00
IKA3160100	R 1826 Plate stand	EA	4	\$98.80
IKA3319000	MS 3 digital vortex mixer 0-3000 rpm	EA	5	\$608.30
IKA3370100	S 10 N - 10 G Dispersing element	EA	4	\$422.75
IKA3386000	R 104 Stand for T10	EA	2	\$70.85

Waterbaths and accessories

Item Code	Description	Unit	Stock on Hand	Price
ISO13-874-29	ISOTEMP 3016S standard ref waterbath 6 L	EA	1	\$3,060.00
ISO13-874-41	ISOTEMP 3013HS heating standard waterbath 13 L	EA	2	\$2,630.00
HAK156-3211	Sahara heated waterbath AC200-S21P 21 L	EA	1	\$2,175.00
HAK170-3050	Waterbath tank plastic type S5P 5 L	EA	1	\$540.00
BAR3161593	Concentric ring covers	EA	1	\$650.00
BAR3166200	Flask clip 1000 mL for precision shaking baths	EA	4	\$142.00

Acids

Item Code	Description	Unit	Stock on Hand	Price
BDH20252.420	Hydrochloric acid 37% NORMAPUR AR	2.5 L	5	\$24.10
FSBH/1111/PB17	Hydrochloric acid MOS 36%	2.5 L	3	\$25.00
FSBP/1280/PB17	Perchloric acid 70% CertiFied AR	2.5 L	36	\$151.90
FSBS/9200/PB17	Sulphuric acid 95-98% low in nitrogen	2.5 L	36	\$25.40
AJA2-500ML	Acetic acid glacial UNILAB	500 mL	5	\$16.00

Bioreagents

Item Code	Description	Unit	Stock on Hand	Price
APPA0839.0010	Ampicillin sodium salt	10 g	1	\$55.30
APPA3261.0100	Dextran sulfate 40 sodium salt BioChemica	100 g	5	\$432.70
FSBBP229-1	Glycerol molecular biology	1 L	23	\$37.10
FSBBP229-4	Glycerol molecular biology	4 L	3	\$114.40
APPA1112.0100	SDS Ultrapure	100 g	18	\$25.60
FSBBP150-20	TEMED Electrophoresis grade	20 g	2	\$18.40
FSBBP154-1	Tris base triple crystallized	1 kg	4	\$102.40
FSBBP169-212	Urea for molecular biology	2.5 kg	3	\$93.70
FSBBP561-1	Water RNA grade sterile	1 L	21	\$31.90
APPA1552.0050	Yeast extract BioChemica	50 g	17	\$10.00

Dries

Item Code	Description	Unit	Stock on Hand	Price
AJA1692-1KG	Aluminium hydroxide UNILAB	1 kg	2	\$71.30
BDH100103M	Aluminium sulphate 18-hydrate AnalaR	500 g	3	\$84.90
BSPAF131.500	Ammonium iron (II) sulphate Pronalys AR	500 g	2	\$22.00
AJA39-5KG	Ammonium iron (II) sulphate UNIVAR	5 kg	1	\$200.90
AJA80-500G	Barium acetate UNIVAR	500 g	1	\$27.80
ACR215410010	Benzalkonium chloride alkyl	1 kg	2	\$175.00
AJA907-500G	Boric acid LABCHEM	500 g	3	\$15.20
AJA102-25KG	Boric acid powder UNILAB	25 kg	6	\$243.70
CMSC018-500G	Calcium sulphate hemihydrate TG	500 g	2	\$7.30
ACR165850250	Cholesteryl oleyl carbonate 99+%	25 g	2	\$52.30
AJA988-100G	Cobalt (II) sulphate UNIVAR	100 g	2	\$145.70
BSPCV180.025	Crystal Violet staining powder	25 g	3	\$8.50
ACR257890010	D(-)-Tagatose 99+%	1 g	2	\$127.20
AJA1364-500G	D-Glucose anhydrous UNILAB	500 g	20	\$9.60
FSBA/5360/53	di-Ammonium hydrogen orthophosphate AR	500 g	33	\$73.10
FSBP/5248/53	di-Potassium hydrogen orthophosphate LR	500 g	45	\$45.70

AJA1234-5KG	di-Sodium hydrogen orthophosphate anhydrous	5 kg	1	\$198.40
AJA358-15KG	di-Sodium hydrogen orthophosphate 12H ₂ O UNILAB	15 kg	2	\$230.00
FSBD/0450/53	EDTA Specified LR	500 g	25	\$30.30
FSBE/0300/53	Eschkas Mixture Certified AR	500 g	3	\$294.40
BSPEL875.20	Ethyl acetate Pronalys AR ACS	20 L	3	\$88.30
AJA1798-1KG	Haematoxylin LABCHEM	1 kg	4	\$866.00
FSBI/0450/53	Iodine Specified LR	500 g	3	\$102.00
AJA1539-500G	Lanthanum oxide UNILAB	500 g	2	\$110.50
CMSAP022-5KG	L-Ascorbic acid BP	5 kg	1	\$20.00
AJA1091-25G	L-Histidine monohydrochloride LABCHEM	25 g	3	\$21.90
AJA624-500G	Lithium sulphate UNIVAR	500 g	6	\$37.70
APPA3445.0025	L-Tryptophan	25 g	2	\$58.50
ACR140641000	L-Tyrosine 99+%	100 g	1	\$73.10
AJA297-5KG	Magnesium chloride UNILAB	5 kg	2	\$45.80
AJA296-5KG	Magnesium chloride UNIVAR	5 kg	1	\$59.00
AJA828-100G	Magnesium oxide UNIVAR	100 g	3	\$15.50
AJA793-500G	Manganese dioxide precipitated UNILAB	500 g	63	\$11.40
AJA598-500G	Manganese dioxide TECHNICAL	500 g	16	\$10.50
AJA1129-100G	Mercury (II) acetate UNILAB	100 g	7	\$36.90
AJA314-100G	Mercury (II) iodide red UNILAB	100 g	1	\$20.30
FSBN/1650/50	Nickel (II) chloride hexahydrate certified AR	250 g	2	\$70.30
ACR388450250	Phosphorus oxybromide 95%	25 g	4	\$118.90
FSBP/4280/65	Potassium chloride certified AR	5 kg	3	\$198.00
AJA384-25KG	Potassium chloride UNILAB	25 kg	4	\$212.50
FSBP/4800/60	Potassium dihydrogen orthophosphate AR	1 kg	3	\$27.50
AJA1079-10KG	Potassium fluoride UNILAB	10 kg	2	\$214.70
FSBP/4880/60	Potassium hexacyanoferrate (III) certified AR	1 kg	7	\$111.50
AJA1202-100G	Silver sulphate UNILAB	100 G	4	\$154.20
AJA680-25KG	Sodium acetate anhydrous UNILAB	25 kg	5	\$186.00
AJA1225-1KG	Sodium carbonate decahydrate UNILAB	1 kg	1	\$31.90
AJA834-500G	Sodium cobaltinitrite UNIVAR	500 g	24	\$81.40
FSBS/3760/60	Sodium dihydrogen orthophosphate 2H ₂ O AR	1 kg	41	\$33.30
AJA471-5KG	Sodium dihydrogen orthophosphate 2H ₂ O UNIVAR	5 kg	41	\$72.90
ACR149190010	Sodium formaldehydesulphoxylate dihydrate 98%	1 kg	2	\$57.90
FSBS/4240/53	Sodium hydrogen carbonate Certified AR	500 g	4	\$12.80
BDH27955.238	Sodium nitrate NORMAPUR AR	250 g	5	\$22.80
AJA2387-25G	Sodium pyruvate LABCHEM	25 g	13	\$24.90

APPA3537.2500	Sodium silicate Solution Pure	2.5 L	1	\$40.20
FSBS/6640/60	Sodium sulphate anhydrous certified AR	1 kg	2	\$36.70
CMSSL066-500G	Sodium sulphate decahydrate LR	500 g	2	\$17.30
BDH102734W	Succinic acid AnalaR	500 g	28	\$65.10
AJA3263-25G	Sudan III C.I. 26100 LABCHEM	25 g	1	\$22.60
AJA1284-500G	Talc TECHNICAL	500 g	22	\$9.30
ACR174110010	Tetraethylene glycol dimethyl ether 99%	1 kg	6	\$99.10
ACR424531000	Thioacetamide 99+% AR ACS	100 g	2	\$194.10
AJA2343-5G	Thymolphthalein LABCHEM	5 g	16	\$10.40
ACR140231000	Trimethyl phosphonoacetate 98%	100 g	3	\$174.50
AJA2311-100G	Tris(hydroxymethyl)methylamine UNIVAR	100 g	2	\$57.60
AJA572-25KG	Urea UNILAB	25 kg	2	\$117.90
ACR116250050	2,2-dimethylglutaric anhydride 98%	5 g	1	\$50.00
ACR430941000	3-aminopropyltriethoxysilane 99% Acroseal	100 mL	2	\$70.90
APPA1421.0001	AEBSF hydrochloride	1 g	3	\$159.30
AJA2297-5G	O-Cresol red LABCHEM	5 g	1	\$52.80

High purity solvents

Item Code	Description	Unit	Stock on Hand	Price
FSBA/0600/PB17	Acetone Certified AR	2.5 L	5	\$14.80
MKT2435-10	Acetone ChromAR for HPLC	4 L	3	\$28.10
FSHA929-1	Acetone OPTIMA grade ACS	1 L	2	\$57.30
BAK9254-03	Acetone ULTRA RESI-ANALYZED	4 L	4	\$56.50
BAK9829-03	Acetonitrile BAKER ANALYZED LC/MS reagent	4 L	32	\$64.00
FSBA/0636/17	Acetonitrile HPLC far UV	2.5 L	6	\$77.30
FSBA/0630/PB17	Acetonitrile HPLC fluorescence grade	2.5 L	32	\$72.10
FSBC603-4	Chloroform Pesticide GC grade ACS	4 L	48	\$68.80
MKT4879-10	Dichloromethane ChromAR for HPLC	4 L	17	\$55.40
LCNPC1040A-2.5L	Dichloromethane Pesticide grade	2.5 L	201	\$25.50
MKTH485-10	Dichloromethane UltimAR	4 L	12	\$74.30
FSBD159-4	Dimethyl sulphoxide HPLC grade	4 L	29	\$105.80
BAK9304-03	Hexanes (95% n-hexane) HPLC	4 L	10	\$98.00
FSBH/0406/17	Hexanes HPLC grade	2.5 L	33	\$23.80
BAK9830-03	Methanol BAKER ANALYSED LC-MS reagent	4 L	12	\$32.60
FSBA452SK-4	Methanol HPLC Certified ACS safe-cote	4 L	68	\$30.00
BAK9831-03	Water Baker Analysed LC - MS Reagent	4 L	16	\$18.60
FSBW5-1	Water HPLC grade	1 L	6	\$17.50
MKT4443-10	Chloroform ChromAR for HPLC	4 L	9	\$59.10

Solvents

Item Code	Description	Unit	Stock on Hand	Price
AJA747-2.5LGL	Dimethyl sulphoxide UNILAB	2.5 L	2	\$19.20
FSBF/1550/PB17	Formamide Specified LR	2.5 L	10	\$79.60
AJA702-20L	Formamide UNILAB	20 L	3	\$216.20
AJA243-2.5LGL	Glycerol UNILAB	2.5 L	4	\$20.10
FSBA/0560/17	Acetone Specified LR	2.5 L	2	\$16.60
AJA97-2.5LGL	Benzyl alcohol UNILAB BP	2.5 L	7	\$34.20
ACR158922500	Glycerol pure 99+%	250 mL	2	\$18.10
BDH10158BG	Methanol AnalaR PE bottle	2.5 L	244	\$10.50
FSBM/4020/17	Methanol Primar	2.5 L	78	\$33.50
BSPML387.2.5-PL	Methanol R Grade	2.5 L	2	\$12.30
ACR326871000	N,N-Dimethylformamide extra dry 99.8%	100 mL	6	\$62.30
FSBH/0421/PB17	N-Hexane Certified AR Grad	2.5 L	7	\$178.90
AJA712-500ML	Nitromethane UNILAB	500 mL	1	\$58.40
FSBP/1440/17	Petroleum spirit 40-60°C speciFied LR	2.5 L	6	\$10.30
BSPPL865.2.5-PL	Propan-2-ol Pronalys AR ACS	2.5 L	6	\$14.70

Miscellaneous chemicals

Item Code	Description	Unit	Stock on Hand	Price
SPX8260A-I	Alternate internal Standard	1 mL	2	\$36.40
SPX8260A-S	Alternate surrogate standard	1 mL	2	\$36.40
SPXBTEX	BTEX standard	1 mL	2	\$23.20
AJA8010070-50	Buffer tablets pH 7.0 makes 100 mL LABCHEM	50 pcs	2	\$23.80
SPX5511-A	Chlorinated disinfection by-products mix	1 mL	2	\$110.70
ROA0442-200/500ML	Conductivity standard 200 μ S/cm	500 mL	6	\$12.50
ROA0442-58MS/500ML	Conductivity standard 58 mS/cm	500 mL	1	\$12.50
FSBK303-1000	CT-AUTO Kjeldahl tablets	1000 pcs	7	\$242.80
ROA0693/250ML	Gold standard 1000 mg/L	250 mL	1	\$38.90
SPX5511-I	Internal standard bromofluorobenzene 10000 μ g/mL	1 mL	2	\$94.20
ROA2693/1L	Mixed solvents Phenol:Chloroform:isoamyl alcohol 30:19:1	1 L	2	\$202.80
BSPNS/40MG/1L	Nitrogen standard 40 mg/L	1 L	4	\$22.30
M-N90212	Paper indicator pH 8.0-10.0 7 mm x 5 m reel	5 reels	8	\$9.20
AJA7872-100	pH sticks 7-14	100 pcs	697	\$8.90
ROA0732/250ML	Potassium standard 10000 mg/L	250 mL	1	\$52.80
M-N91304	Quantofix Copper test sticks 6 x 95 mm	100 pcs	4	\$26.00
ROA0741/250ML	Scandium standard 1000 mg/L	250 mL	1	\$77.80
ROA1037/ 30BRX/100ML	Sucrose standard 30 brix	100 mL	2	\$17.40

AJA2574-100	Universal indicator cards	100 pcs	4	\$141.20
SPXBTEX-1	VOA Standard matrix methanol-P&T	1 mL	2	\$29.80
SPX8260-BIG-MIX	Volatile organics mix (76 Components)	1 mL	2	\$266.00
SPXS-550	4-bromofluorobenzene 1000 µg/mL	1 mL	2	\$55.30
BSPCOD/1500MG/1L	COD standard 1500 mg/L	1 L	5	\$24.70
ROA0647/1L	Cerium sulphate 0.1N volumetric solution	1 L	1	\$64.00
FSBJ/4350C/90	Hydrochloric acid 0.1 mol/L concentrate	6 pcs	4	\$46.70
AJA2653-20L	Neutral buffered formalin LABCHEM	20 L	29	\$20.50
ROA1183/500ML	Sodium hydroxide 0.1%	500 mL	6	\$16.70
FSBS/4930/25	Sodium hydroxide 46%	25 L	44	\$119.90
BDH180465G	Sulphuric acid 0.05 mol/L 0.1N ConvoL	6 pcs	2	\$79.20
AJA909-2.5LPL	Calcium chloride 40% solution UNIVAR PL	2.5 L	1	\$69.00
FSBC/2120/60	Calcium oxide powder laboratory reagent	1 kg	4	\$18.90
RBSRBS50PF	Multipurpose alkaline cleaning agent	5 L	3	\$85.10
RBSRBST105-5L	RBST105 concentrated liquid detergent	5 L	2	\$72.00
BSPA97.2.5	Buffer pH 4.0 colour coded red	2.5 L	2	\$13.90
ACR351450100	Methyl sulphoxide-d6 with 0.03% TMS	10 mL	2	\$30.30



Better together

Alfa Aesar and Acros Organics, a better source for organics

Find out more thermofisher.com.au/alfaesar

ThermoFisher
SCIENTIFIC

© 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.