



Innovation in tissue processing

unmatched tissue quality

The Thermo Scientific™ Excelsior AS™ Tissue Processor provides exceptional tissue quality with minimal user interaction. You get the superb tissue quality and proven reliability you expect, plus a safer, more efficient working environment.

Enhanced tissue quality

Cassette baskets are gently rotated inside the circular chamber, providing more effective agitation than an impeller.

Unmatched user safety

When the chamber is opened, a fan pulls fumes through a filter to protect users from formaldehyde exposure.

Intelligent reagent management for lowest cost of ownership

Alcohol quality measurement enables you to extend your reagent life and provides significant cost savings.

Streamline workflow and increase productivity

Reduce maintenance time by replacing one bottle of each reagent at a time.



TRADE-IN OFFER

\$55,000

Valid only for trade-in of any tissue processor until December 31st, 2017.

ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-292
To fax an order, use 1800-067-639
Visit us online: www.thermofisher.com.au

©2017 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917