



thermo scientific

Thermo Scientific **Automation Tips**

Quick Reference Guide

ThermoFisher
SCIENTIFIC

Automation tips for a variety of automated liquid handlers

Easily identify your workstation from the chart below to find a complete list of compatible Thermo Scientific automation tips for that instrument:

Table of Contents

Tecan™

Freedom EVO™	2
Genesis™, MiniPrep™, Cavro™	3
Aquarius™, Te-MO™	4

Beckman™

Biomek® FX, Biomek NX.....	5
Multimek™ AP96, Biomek 3000/4000	6
Multimek AP384, Biomek 1000, Biomek 2000, Multimek II	7

Perkin Elmer/Packard™

JANUS™	8
MultiPROBE™ II HT, MultiPROBE II HT EX.....	9
PlateTrak™, MiniTrak™, Evolution P-3™	10

Qiagen™

BioRobot™ 3000/8000/9600/9604, Rosys, Colibri™	11
---	----

Agilent™

VPrep Workstation™, Bravo™ Workstation	11
--	----

Eppendorf™

epMotion™	11
-----------------	----

Caliper/Zymark™

Caliper/Zymark™ RapidPlate™, Zymark Presto™, Zymark SciClone™, Zymark Allegro™, Zymark Zephyr™	12
--	----

BioTek™

Precision™, Precision XS	12
--------------------------------	----

Molecular Devices™

FLIPR™	12
--------------	----

Thermo Scientific™ D.A.R.T.'s™

Versette™ Automated Liquid Handler	13
Matrix™ PlateMate™, Matrix Hdyra.....	14
Matrix SerialMate™	15



Tecan

Freedom EVO

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
LiHa	901-011	10 µL	Conductive	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-021	10 µL	Clear	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-251	20 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-252	20 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-253	20 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	901-261	20 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-262	20 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-011	50 µL	Conductive	Sterile, Filter	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-021	50 µL	Clear	Sterile, Filter	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-251	50 µL	Conductive	Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-252	50 µL	Conductive	Non-Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-253	50 µL	Conductive	Non-Sterile, Econopak	Purple	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	902-261	50 µL	Clear	Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-262	50 µL	Clear	Non-Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-011	175 µL	Conductive	Sterile, Filter	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-011G	175 µL	Conductive	Sterile, Filter, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-021	175 µL	Clear	Sterile, Filter	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-033	175 µL	Conductive	Non-Sterile, Filter, Econopak	Red	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	903-251	200 µL	Conductive	Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-251G	200 µL	Conductive	Sterile, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-252	200 µL	Conductive	Non-Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-252G	200 µL	Conductive	Non-Sterile, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-253	200 µL	Conductive	Non-Sterile, Econopak	Red	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	903-261	200 µL	Clear	Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-262	200 µL	Clear	Non-Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-011	1000 µL	Conductive	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-011G	1000 µL	Conductive	Sterile, Filter, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-013	1000 µL	Conductive	Non-Sterile, Filter, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	904-021	1000 µL	Clear	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-251	1000 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-251G	1000 µL	Conductive	Sterile, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-252	1000 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-252G	1000 µL	Conductive	Non-Sterile, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
904-253	1000 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case	
904-261	1000 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
904-262	1000 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
Te-MO 96	906-021	50 µL	Clear	Sterile, Filter	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	906-261	50 µL	Clear	Sterile	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	906-262	50 µL	Clear	Non-Sterile	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	907-021	100 µL	Clear	Sterile, Filter	Red	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	907-261	100 µL	Clear	Sterile	Red	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	907-262	100 µL	Clear	Non-Sterile	Red	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	908-021	150 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	908-261	200 µL	Clear	Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	908-262	200 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Tecan

Genesis, MiniPrep, Cavro

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
LiHa	901-011	10 µL	Conductive	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-021	10 µL	Clear	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-251	20 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-252	20 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-253	20 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	901-261	20 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-262	20 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-011	50 µL	Conductive	Sterile, Filter	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-021	50 µL	Clear	Sterile, Filter	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-251	50 µL	Conductive	Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-252	50 µL	Conductive	Non-Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-253	50 µL	Conductive	Non-Sterile, Econopak	Purple	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	902-261	50 µL	Clear	Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-262	50 µL	Clear	Non-Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-011	175 µL	Conductive	Sterile, Filter	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-011G	175 µL	Conductive	Sterile, Filter, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-021	175 µL	Clear	Sterile, Filter	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-033	175 µL	Conductive	Non-Sterile, Filter, Econopak	Red	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	903-251	200 µL	Conductive	Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-251G	200 µL	Conductive	Sterile, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-252	200 µL	Conductive	Non-Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-252G	200 µL	Conductive	Non-Sterile, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-253	200 µL	Conductive	Non-Sterile, Econopak	Red	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	903-261	200 µL	Clear	Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-262	200 µL	Clear	Non-Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-011	1000 µL	Conductive	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-011G	1000 µL	Conductive	Sterile, Filter, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-013	1000 µL	Conductive	Non-Sterile, Filter, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	904-021	1000 µL	Clear	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-251	1000 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-251G	1000 µL	Conductive	Sterile, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-252	1000 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-252G	1000 µL	Conductive	Non-Sterile, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-253	1000 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	904-261	1000 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-262	1000 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack

Tecan

Aquarius, Te-MO

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
Te-MO 96	906-021	50 µL	Clear	Sterile, Filter	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	906-261	50 µL	Clear	Sterile	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	906-262	50 µL	Clear	Non-Sterile	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	907-021	100 µL	Clear	Sterile, Filter	Red	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	907-261	100 µL	Clear	Sterile	Red	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	907-262	100 µL	Clear	Non-Sterile	Red	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	908-021	150 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	908-261	200 µL	Clear	Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	908-262	200 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Beckman

Biomek FX, Biomek NX

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
AP96	912-021	10 µL	Clear	Sterile, Filter	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	912-261	10 µL	Clear	Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	912-262	10 µL	Clear	Non-Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	916-021	40 µL	Clear	Sterile, Filter	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	916-261	50 µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	916-262	50 µL	Clear	Non-Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	916-262-05	50 µL	Clear	Non-Sterile, Low Retention	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-021	20 µL	Clear	Sterile, Filter	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-021-05	20 µL	Clear	Sterile, Filter, Low Retention	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-261	20 µL	Clear	Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-261-05	20 µL	Clear	Sterile, Low Retention	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-262	20 µL	Clear	Non-Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-262-05	20 µL	Clear	Non-Sterile, Low Retention	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-021	130 µL	Clear	Sterile, Filter	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-021-05	130 µL	Clear	Sterile, Filter, Low Retention	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-021G	130 µL	Clear	Sterile, Filter, Wide Bore	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-261	250 µL	Clear	Sterile	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-261-05	250 µL	Clear	Sterile, Low Retention	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-261G	250 µL	Clear	Sterile, Wide Bore	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-262	250 µL	Clear	Non-Sterile	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
919-262-05	250 µL	Clear	Non-Sterile, Low Retention	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
919-262G	250 µL	Clear	Non-Sterile, Wide Bore	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
AP384	935-261	30 uL	Clear	Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	935-262	30 uL	Clear	Non-Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Beckman

Multimek AP96, Beckman Biomek 3000/4000

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
AP96	912-021	10 µL	Clear	Sterile, Filter	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	912-261	10 µL	Clear	Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	912-262	10 µL	Clear	Non-Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	916-021	40 µL	Clear	Sterile, Filter	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	916-261	50 µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	916-262	50 µL	Clear	Non-Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	916-262-05	50 µL	Clear	Non-Sterile, Low Retention	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-021	20 µL	Clear	Sterile, Filter	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-021-05	20 µL	Clear	Sterile, Filter, Low Retention	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-261	20 µL	Clear	Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-261-05	20 µL	Clear	Sterile, Low Retention	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-262	20 µL	Clear	Non-Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	918-262-05	20 µL	Clear	Non-Sterile, Low Retention	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-021	130 µL	Clear	Sterile, Filter	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-021-05	130 µL	Clear	Sterile, Filter, Low Retention	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-021G	130 µL	Clear	Sterile, Filter, Wide Bore	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-261	250 µL	Clear	Sterile	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-261-05	250 µL	Clear	Sterile, Low Retention	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-261G	250 µL	Clear	Sterile, Wide Bore	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-262	250 µL	Clear	Non-Sterile	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-262-05	250 µL	Clear	Non-Sterile, Low Retention	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	919-262G	250 µL	Clear	Non-Sterile, Wide Bore	Light Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Beckman

Multimek AP384

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
AP384	935-261	30 µL	Clear	Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	935-262	30 µL	Clear	Non-Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Beckman

Biomek 1000, Biomek 2000

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
MP20	915-021	20 µL	Clear	Sterile, Filter	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	915-261	20 µL	Clear	Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	915-262	20 µL	Clear	Non-Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
MP200	917-021	130 µL	Clear	Sterile, Filter	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	917-261	200 µL	Clear	Sterile	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	917-262	200 µL	Clear	Non-Sterile	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Beckman

Multimek II

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
Type II	2760	50 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	2790	175 µL	Clear	Sterile, Filter	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	3762	50 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	3771-05	300 µL	Clear	Sterile, Low Retention	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	3772-05	300 µL	Clear	Non-Sterile, Low Retention	Green	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Perkin Elmer/Packard

JANUS

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging	
Varispan	901-251	20 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	901-252	20 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	901-253	20 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case	
	901-261	20 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	901-262	20 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	902-251	50 µL	Conductive	Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	902-252	50 µL	Conductive	Non-Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	902-253	50 µL	Conductive	Non-Sterile, Econopak	Purple	Hanging Trays	96 Tips/Tray, 50 Trays/Case	
	902-261	50 µL	Clear	Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	902-262	50 µL	Clear	Non-Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	903-011	175 µL	Conductive	Sterile, Filter	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	903-021	175 µL	Clear	Sterile, Filter	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	903-251	200 µL	Conductive	Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	903-251G	200 µL	Conductive	Sterile, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	903-252	200 µL	Conductive	Non-Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	903-252G	200 µL	Conductive	Non-Sterile, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	903-253	200 µL	Conductive	Non-Sterile, Econopak	Red	Hanging Trays	96 Tips/Tray, 50 Trays/Case	
	903-261	200 µL	Clear	Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	903-262	200 µL	Clear	Non-Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	904-251	1000 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	904-251G	1000 µL	Conductive	Sterile, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	904-252	1000 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	904-252G	1000 µL	Conductive	Non-Sterile, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	904-253	1000 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case	
	904-261	1000 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	904-262	1000 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	904-261	1000 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	904-262	1000 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack	
	MDT 96	925-021	20 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
		925-261	20 µL	Clear	Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
925-262		20 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
927-021		50 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
927-261		50 µL	Clear	Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
927-261-05		50 µL	Clear	Sterile, Low Retention	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
927-262		50 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
929-021		235 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
929-261		235 µL	Clear	Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
929-262		235 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
MDT 384	931-021	25 µL	Clear	Sterile, Filter	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
	931-261	30 µL	Clear	Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	
	931-262	30 µL	Clear	Non-Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case	

Perkin Elmer/Packard

MultiPROBE II HT, MultiPROBE II HT EX

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
Varispan	901-251	20 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-252	20 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-253	20 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	901-261	20 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	901-262	20 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-251	50 µL	Conductive	Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-252	50 µL	Conductive	Non-Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-253	50 µL	Conductive	Non-Sterile, Econopak	Purple	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	902-261	50 µL	Clear	Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	902-262	50 µL	Clear	Non-Sterile	Purple	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-011	175 µL	Conductive	Sterile, Filter	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-021	175 µL	Clear	Sterile, Filter	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-251	200 µL	Conductive	Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-251G	200 µL	Conductive	Sterile, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-252	200 µL	Conductive	Non-Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-252G	200 µL	Conductive	Non-Sterile, Wide Bore	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-253	200 µL	Conductive	Non-Sterile, Econopak	Red	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	903-261	200 µL	Clear	Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	903-262	200 µL	Clear	Non-Sterile	Red	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-251	1000 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-251G	1000 µL	Conductive	Sterile, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-252	1000 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-252G	1000 µL	Conductive	Non-Sterile, Wide Bore	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-253	1000 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	904-261	1000 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-262	1000 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-261	1000 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	904-262	1000 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack

Perkin Elmer/Packard

PlateTrak, MiniTrak, Evolution P-3

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
MDT 96	925-021	20 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	925-261	20 µL	Clear	Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	925-262	20 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	927-021	50 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	927-261	50 µL	Clear	Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	927-261-05	50 µL	Clear	Non-Sterile, Low Retention	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	927-262	50 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	929-021	235 µL	Clear	Sterile, Filter	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	929-261	235 µL	Clear	Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	929-262	235 µL	Clear	Non-Sterile	Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
MDT 384	931-021	25 µL	Clear	Sterile, Filter	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	931-261	30 µL	Clear	Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	931-262	30 µL	Clear	Non-Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Qiagen

BioRobot 3000/8000/9600/9604, Rosys, Colibri

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
Independent Channel	951-011	250 µL	Conductive	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	951-021	250 µL	Clear	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	951-251	300 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	951-252	300 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	951-253	300 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	951-261	300 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	951-262	300 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	953-011	1000 µL	Conductive	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	953-021	1000 µL	Clear	Sterile, Filter	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	953-251	1100 µL	Conductive	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	953-252	1100 µL	Conductive	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	953-253	1100 µL	Conductive	Non-Sterile, Econopak	Gray	Hanging Trays	96 Tips/Tray, 50 Trays/Case
	953-261	1100 µL	Clear	Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack
	953-262	1100 µL	Clear	Non-Sterile	Gray	Hanging Trays	96 Tips/Tray, 10 Trays/Pack

Agilent

VPrep Workstation, Bravo Workstation

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
96ST, 384ST	946-261	10 µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	946-262	10 µL	Clear	Non-Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	947-261	30 µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	947-262	30 µL	Clear	Non-Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	948-261	50 µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	948-262	50 µL	Clear	Non-Sterile	Black	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Eppendorf

epMotion Workstations

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
Independent Channel	965-021	50 µL	Clear	Sterile, Filter	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case
	965-261	50 µL	Clear	Sterile	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case
	965-262	50 µL	Clear	Non-Sterile	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case
	966-021	300 µL	Clear	Sterile, Filter	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case
	966-261	300 µL	Clear	Sterile	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case
	966-262	300 µL	Clear	Non-Sterile	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case
	967-021	1000 µL	Clear	Sterile, Filter	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case
	967-261	1000 µL	Clear	Sterile	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case
	967-262	1000 µL	Clear	Non-Sterile	White	SBS	96 Tips/Rack, 10 Racks/Pack, 4 Packs/Case

Caliper/Zymark

Caliper/Zymark RapidPlate, Caliper Presto, Caliper SciClone, Caliper Allegro, Caliper Zephyr

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
96-well	920-021	50 µL	Clear	Sterile, Filter	Green	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	920-261	50 µL	Clear	Sterile	Green	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	920-262	50 µL	Clear	Non-Sterile	Green	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	921-021	20 µL	Clear	Sterile, Filter	Green	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	921-261	100 µL	Clear	Sterile	Green	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	921-262	100 µL	Clear	Non-Sterile	Green	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	923-021	100 µL	Clear	Sterile, Filter	Rose	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	923-261	200 µL	Clear	Sterile	Rose	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	923-262	200 µL	Clear	Non-Sterile	Rose	Non-SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

BioTek

Precision, Precision XS

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
Independent Channel	984-021	50 µL	Clear	Sterile, Filter	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	984-261	50 µL	Clear	Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	984-262	50 µL	Clear	Non-Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	985-021	80 µL	Clear	Sterile, Filter	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	985-261	100 µL	Clear	Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	985-262	100 µL	Clear	Non-Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	986-021	150 µL	Clear	Sterile, Filter	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	986-261	200 µL	Clear	Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	986-262	200 µL	Clear	Non-Sterile	Light Blue	SBS	96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Molecular Devices

FLIPR

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
Type A CCS	931-021	25 µL	Clear	Sterile, Filter	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	931-261	30 µL	Clear	Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	931-262	30 µL	Clear	Non-Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
Type B Liberty	937-261	30 µL	Clear	Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case
	937-262	30 µL	Clear	Non-Sterile	Black	SBS	384 Tips/Rack, 10 Racks/Pack, 5 Packs/Case

Thermo Scientific D.A.R.T.'s

Versette Automated Liquid Handler

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
96-well	5501	25µL	Clear	Non-Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5502-11	25µL	Clear	Sterile	Black	Magazines	96 Tips/Mgz, 960 Tips/Case
	5503-11	30µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5506-11	30µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5507	25µL	Clear	Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5508	30µL	Clear	Sterile, Filter, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5516-11	300µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5517-11	300µL	Clear	Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5518-11	300µL	Clear	Sterile, Filter, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5526	100µL	Clear	Non-Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5527	100µL	Clear	Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5528	100µL	Clear	Sterile, Filter	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5536	300µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5537	300µL	Clear	Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5538	300µL	Clear	Sterile, Filter, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5546	300µL	Clear	Non-Sterile, Wide Bore, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5547	300µL	Clear	Sterile, Wide Bore, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5548	300µL	Clear	Sterile, Filter, Wide Bore, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5562-11	220µL	Clear	Sterile, Extended Length	Black	SBS	96 Tips/Rack, 960 Tips/Case
	5586	30µL	Clear	Non-Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5587	30µL	Clear	Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5588	30µL	Clear	Sterile, Filter	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5606-11	12.5/30µL	Clear	Non-Sterile	Black	SBS	96 Tips/Rack, 960 Tips/Case
	5607	12.5/30µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 960 Tips/Case
	5608	12.5/30µL	Clear	Sterile, Filter	Black	SBS	96 Tips/Rack, 960 Tips/Case
	5618	1205µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 960 Tips/Case
	5706	125/300µL	Clear	Non-Sterile	Black	SBS	96 Tips/Rack, 960 Tips/Case
	5707	125/300µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 960 Tips/Case
	5708	125/300µL	Clear	Sterile, Filter	Black	SBS	96 Tips/Rack, 960 Tips/Case
	5806	1250µL	Clear	Non-Sterile	Black	SBS	96 Tips/Rack, 768 Tips/Case
	5807	1250µL	Clear	Sterile	Black	SBS	96 Tips/Rack, 768 Tips/Case
	5808	1250µL	Clear	Sterile, Filter	Black	SBS	96 Tips/Rack, 768 Tips/Case
	384-well	5301	12.5µL	Clear	Non-Sterile	Black	SBS
5302		12.5µL	Clear	Sterile	Black	SBS	384 Tips/Rack, 3840 Tips/Case
5316		30µL	Clear	Non-Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5317		30µL	Clear	Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5318		30µL	Clear	Sterile, Filter	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5321		100µL	Clear	Non-Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5326		100µL	Clear	Non-Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5327		100µL	Clear	Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5328		100µL	Clear	Sterile, Filter	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5416		30µL	Clear	Non-Sterile, Extended Length	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5416R		30µL	Clear	Non-Sterile, Extended Length	Black	SBS	384 Tips/Rack, 7680 Tips/Case
5417		30µL	Clear	Sterile, Extended Length	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
5417R		30µL	Clear	Sterile, Extended Length	Black	SBS	384 Tips/Rack, 7680 Tips/Case
5418		30µL	Clear	Sterile, Filter, Extended Length	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case

Thermo Scientific D.A.R.T.'s

Matrix PlateMate, Matrix Hydra

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
96-well	5501	25µL	Clear	Non-Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5502-11	25µL	Clear	Sterile	Black	Magazines	96 Tips/Mgz, 960 Tips/Case
	5503-11	30µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5506-11	30µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5507	25µL	Clear	Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5508	30µL	Clear	Sterile, Filter, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5516-11	300µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5517-11	300µL	Clear	Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5518-11	300µL	Clear	Sterile, Filter, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5526	100µL	Clear	Non-Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5527	100µL	Clear	Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5528	100µL	Clear	Sterile, Filter	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5536	300µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5537	300µL	Clear	Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5538	300µL	Clear	Sterile, Filter, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5546	300µL	Clear	Non-Sterile, Wide Bore, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5547	300µL	Clear	Sterile, Wide Bore, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5548	300µL	Clear	Sterile, Filter, Wide Bore, Extended Length	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5586	30µL	Clear	Non-Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5587	30µL	Clear	Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
5588	30µL	Clear	Sterile, Filter	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case	
	5562-11	220µL	Clear	Sterile, Extended Length	Black	SBS	96 Tips/Rack, 960 Tips/Case
384-well	5301	12.5µL	Clear	Non-Sterile	Black	SBS	384 Tips/Rack, 3840 Tips/Case
	5302	12.5µL	Clear	Sterile	Black	SBS	384 Tips/Rack, 3840 Tips/Case
	5316	30µL	Clear	Non-Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5317	30µL	Clear	Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5318	30µL	Clear	Sterile, Filter	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5321	100µL	Clear	Non-Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5326	100µL	Clear	Non-Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5327	100µL	Clear	Sterile	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5328	100µL	Clear	Sterile, Filter	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5416	30µL	Clear	Non-Sterile, Extended Length	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5416R	30µL	Clear	Non-Sterile, Extended Length	Black	SBS	384 Tips/Rack, 7680 Tips/Case
	5417	30µL	Clear	Sterile, Extended Length	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case
	5417R	30µL	Clear	Sterile, Extended Length	Black	SBS	384 Tips/Rack, 7680 Tips/Case
	5418	30µL	Clear	Sterile, Filter, Extended Length	Black	Magazines	384 Tips/Mgz, 7680 Tips/Case

Thermo Scientific D.A.R.T.'s

Matrix SerialMate

Head	Thermo Scientific Cat. Number	Volume	Clear/Conductive	Description	Rack Colors	Rack Format	Packaging
96-well	5541	250µL	Clear	Non-Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5542	250µL	Clear	Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5545	250µL	Clear	Sterile, Filter	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5561	220µL	Clear	Non-Sterile, Extended Length	Black	Magazines	96 Tips/Mgz, 960 Tips/Case
	5571	30µL	Clear	Non-Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
	5572	30µL	Clear	Sterile	Black	Magazines	96 Tips/Mgz, 1920 Tips/Case
384-well	5311	30µL	Clear	Non-Sterile	Black	SBS	384 Tips/Rack, 3840 Tips/Case
	5312	30µL	Clear	Sterile	Black	SBS	384 Tips/Rack, 3840 Tips/Case

In Australia:

For customer service, call 1300-735-292

To fax an order, use 1800-067-639

To email an order, ordersau@thermofisher.com

Find out more at

thermofisher.com.au/automation

In New Zealand:

For customer service, call 0800-933-966

To fax an order, use 0800-329-246

To email an order, ordersnz@thermofisher.com

Find out more at

thermofisher.co.nz/automation

ThermoFisher
SCIENTIFIC

© 2019 Thermo Fisher Scientific Inc. All rights reserved. Please consult your local sales representative for details. "Bravo" and "VPrep" (Agilent Technologies); Beckman, "Biomek" and "Multimek" (Beckman Coulter); "BioTek", "Precision" and "Precision XS" (BioTek, Inc.); "Colibri" (Colibri Robotics, Inc.); "Eppendorf" and "epMotion" (Eppendorf AG); "Hamilton" and "Star", (Hamilton Robotics); "FLIPR" and "Molecular Devices" (Molecular Devices); "JANUS", "PerkinElmer", "Packard", "MultiPROBE" and "PlateTrak", "MiniTrak", "Evolution P-3" and "Varispan" (Perkin Elmer); "Tecan", "Genesis", "Freedom Evo", "MiniPrep", "Cavro", "Te-MO" and "Aquarius" (Tecan Group); "BioRobot", "Qiagen", and "Rosys" (Qiagen GmbH); "Allegro", "Caliper", "Presto", "Rapid Plate", "SciClone", "Zephyr" and "Zymark" (Caliper Life Sciences) are trademarks of their respective owners. All other trademarks owned by Thermo Fisher Scientific or its subsidiaries.